

Granada

Model No. **S743**

ROMA



Special Features

- Attractive removeable modety panel
- Full LED lights set
- LCD dashboard for greater clarity
- Dual mode wig-wag lever
- Automatic speed reduction on full lock turns in 8mph setting for safety
- Flip up and width and height adjustable armrests
- Infinite angle tiller adjustment
- Rear disc brake system
- Non marking pneumatic tyres
- Comfortable height adjustable, reclining seat
- Large front basket with carry handle
- Automatic brake lights
- Full adjustable suspension

Specifications



Model No.*	S743	Maximum User Weight	160 kg / 25 st
Class	3	Battery Size 12v x 2	50Ah
Overall Length	1300 mm / 51"	Range**	Up to 40 km / 25 miles
Overall Width	660 mm / 26"	Overall Weight	110 kgs / 242 lbs
Overall Height (folded, no seat)	640 mm / 25"	Maximum Speed	6 & 12 kmph / 4 & 8 mph
Overall Height (top of mirror)	1240 mm / 49"	Wheel / Tyre Size (Front)	300 mm / 12"
Seat Width	460 mm / 18"	(Back)	300 mm / 12"
Seat Depth	410 mm / 16"	Turning Radius	1265 mm / 50"
Seat Height (to ground)	600 mm-720 mm / 24"-28.5"	Battery Charger	5A
Seat Height (to floorpan)	400 mm-520 mm / 16"-20.5"	Ground Clearance	120 mm / 5"
Backrest Height (from seat)	390 mm / 15"		
Safe Working Slope	8° / 14%		

* Lower limb amputees should refer to user manual before use. ** Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel. *** Measurement to underside of chassis.