


Rascal Pioneer



Optional Extra Back Box



OPTIONAL ACCESSORY



Waterproof Panel and Twin Mirrors



Rascal Pioneer

Rascal Pioneer Specification:

Type Class	C / 3
Overall length	147cm / 58"
Overall width	64.5cm / 25"
Floor to seat height	65-68cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	159kgs / 25 stone
Battery voltage and capacity	12V / 75Ah x 2**
Maximum speed	12.9km/h / 8mph
Ground clearance	9cm / 3.5"
Maximum kerb height	10cm / 4"
Range*	up to 56km / 35 miles
The maximum safe slope	10°
Turn-around width	2.02m / 6'7"
Colour	Black / Red

*Range on full charge and flat ground.

**Free upgrade to 80Ah batteries currently available.

Main Features:

- Front and rear coil suspension
- Removable front basket
- Articulated tiller adjustment with useful battery charging port
- Modern delta handlebars
- Waterproof panel and twin mirrors
- Large shock absorbing tyres for maximum comfort on the move
- Sliding swivel seat with armrests and headrest for maximum comfort

Optional Extras:

- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat
- An exclusive Rascal basket bag



Rascal Pioneer Profile:

The Rascal Pioneer is a proven performer with a fantastic reliability and durability record. With full coil suspension and a 35 mile range its powerful motor makes it an ideal scooter for the more challenging environments. It features a delta-bar tiller with hydraulic adjustment, a sliding swivel seat with headrest and padded arms as well as waterproof switches and a large basket making the Pioneer one of the most popular scooters on the market.



Local Rascal Mobility Retailer:

E. and O.E- v 5 Final-July 2015