

Reno II

Model No. **P319**

ROMA



Special Features

- UK designed & registered by Roma.
- New battery pack 2 x 22ah.
- New patented docking system.
- Intelligent, regenerative, electromagnetic brakes.
- New frame locking mechanism.
- Large storage basket.
- On and Off board charger
- Easily disassembly without tools for carrying, storage and transportation.
- Compact and lightweight design with excellent manoeuvrability.
- Simplified frame structure for easy maintenance access.
- Non-marking puncture proof tyres.

Specifications



Model No.*	P319	Battery Size 12v x 2	22 Ah
Class (Invalid Carriage)	2	Range**	20 km / 12.5 miles
Overall Length	869 mm / 35.2"	Motor	DC24V, 90W
Overall Width	571 mm / 22.5"	Maximum Speed	4 mph / 6.4 kmph
Overall Height	860 mm / 34"	Wheel Size (Front)	150 mm / 6"
Overall Height (Seat removed)	427 mm / 16.8"	Wheel Size (Back)	230 mm / 9"
Seat Width	460 mm / 18"	Turning Radius	701 mm / 27.6"
Seat Depth	380 mm / 15"	Battery Charger	4 Amp Off Board
Seat Height (to ground)	532 mm / 21"	Battery Pack Weight	14 Kgs / 31 lbs
Seat Height (to floorpan)	427 mm / 16.8"		
Safe Working Slope	6° / 10%		
Maximum User Weight	136 kgs / 21st		

* Lower limb amputees should refer to user manual before use. ** Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.