

ROMA



Marbella

Model No. **P200**

Special Features

- Adjustable angle seat and back
- Adjustable height and width arms
- Twin detachable battery packs
- Backrest folds down for transport and storage
- Heaviest part only 18 kgs (main chassis)
- Crash tested



Specifications



Model No.*	P200	Charger	5 amp off board
Class	2	Range**	Up to 40 km / 25 miles
Overall Length	1066 mm / 42"	Overall Weight	80 kg / 176 lbs
Overall Width	610 mm / 24"	Maximum Speed	6 kmph / 4 mph
Overall Height	960 mm / 38"	Wheel / Tyre size	360 mm / 14"
Seat Width	460 mm / 18"	Turning Radius	580 mm / 23"
Seat Depth	400 mm / 16"	Ground Clearance	100 mm / 4"
Seat Height	495-520 mm / 19 - 20"	Colour	Red
Backrest Height from Seat	430 mm / 17"	Crash Tested	
Safe Working slope	8° / 12%		
Maximum User Weight	136 kg / 21 st		
Battery Size	2 x 12v 50 amp		

* Lower limb amputees should refer to user manual before use. ** Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.