

SHOPRIDER



# Grande

Model No. **TE-9AS**

## Special Features

- Full LED lights set
- Full suspension
- 50 amp power battery
- High and low speed settings
- Fully adjustable captain seat
- Black non-marking wide profile pneumatic tyres
- Free wheel release handle
- Delta handle
- USB charging port
- Flip up adjustable armrests
- Tiller angle adjustment lever
- Rear indicators & hazard lights



## Specifications



Model No.*	<b>TE-9AS</b>	Safe Working Slope	<b>8° / 14%</b>
Class (Invalid Carriage)	<b>3</b>	Maximum User Weight	<b>159 kgs / 25 st</b>
Overall Length	<b>1400 mm / 55"</b>	Battery Size 12v x 2	<b>36 Ah / 50 Ah (Option)</b>
Overall Width	<b>600 mm / 24"</b>	Range**	<b>28 km / 18 miles</b>
Overall Height (folded, no seat)	<b>1360 mm / 54"</b>	Weight	<b>99 kgs / 218 lbs</b>
Overall Height (top of mirror)	<b>1170 mm / 46"</b>	Maximum Speed	<b>12 mph / 7.5 kmph</b>
Seat Width	<b>460 mm / 18"</b>	Wheel / Tyre Size (Front)	<b>330 mm / 13"</b>
Seat Depth	<b>430 mm / 17"</b>	Wheel / Tyre Size (Back)	<b>330 mm / 13"</b>
Seat Height (to ground)	<b>660 mm / 26"</b>	Turning Radius	<b>1680 mm / 66"</b>
Armrest Height	<b>170 mm / 7"</b>	Battery Charger	<b>5 Amp Standard</b>
Backrest Height (from Seat)	<b>730 mm / 29"</b>	Ground Clearance***	<b>10 mm / 4"</b>
Seat Height (to Floorpan)	<b>470 mm / 18.5"</b>		

\* Lower limb amputees should refer to user manual before use. \*\* Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.