



WHERE YOU WANT TO BE.

Little Star





Little Star

Quality and style
for everyone





Looking for the ideal portable scooter which combines simplicity of use and style?

Try the bright Little Star from Sterling.



Simple controls. Small, easy to use controller with the added security of an off/on key switch




Footmat for comfort



Folds easily for transportation or storage



Lightweight battery pack
7.3kg (16lbs)




Easy access with simple to use flip back armrests. Padded seat and backrest for maximum comfort



Narrow profile with a compact width of 50cm (20") navigating through doorways causes no problem

Little Star



Technical data sheet

Version	4 Wheels	Heaviest part	20kg (44lbs)
Speed	4mph (6kph)	Turning radius	110cm (43.3")
Batteries	10Ah	Maximum safe slope	6°
Maximum range	7.4 miles (12km)	Ground clearance	3.5cm (1.3")
Colour	Black	Front wheel size	18cm (7")
Maximum user weight	113kg (17.8 stone)	Rear wheel size	19cm (7.5")
Total weight	35kg (77lbs)	Lights and indicators package	No

Find out more today!

Visit www.sterlingscooters.co.uk for more information, to find your local dealer or to request a demonstration.

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively is available on request in large text.



Sunrise Medical Limited

Thorns Road
Brierley Hill
West Midlands
DY5 2LD

Tel: +44 (0) 845 605 66 88
Fax: +44 (0) 845 605 66 89

www.SunriseMedical.co.uk

* Range given represents test conditions. Please note that temperature, user mass, terrain and battery capacity may affect the overall performance.