

# TDX<sup>®</sup>2 **ULTRA**

Retail Order Form  
Max user weight 136kg (21.4 stone)

Featuring

**ULTRALow maxx**  
By Motion Concepts

- Class leading shear reduction with our ESR Mechanism
- Achieve the perfect fit with simple and precise adjustment features
- Maximise comfort and pressure relief with the advanced tilt and recline
- Superior drive performance with Stability Lock and SureStep<sup>®</sup> Technology



Invacare Ltd  
 Pencoed Technology Park  
 Pencoed CF35 5AQ  
 Tel: 01656 776222 Fax: 01656 776220  
 Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)

Power Wheelchair  
 Dealer Prescription Form  
 LPFTDX2ULTRA010116DEALER  
 January 2016

E-spares cat Now Available Online  
 at [www.invacare.co.uk](http://www.invacare.co.uk)

Customer Name :  
 Dealer :  
 Address :  
 Tel No:  
 Delivery address :

Account No :  
 Order No :  
 Order Qty :  
 Order Date :

### Standard Features



- Ultra Low Maxx seating by Motion Concepts
- Max user weight 136kg (21.4 stone)
- Compact, manoeuvrable, CWD TDX SP2 chassis
- Seat widths of 405-610mm (16"- 24"), adjustable between 405-505mm (16"- 20") and 505-610mm (20"-24").
- Seat depth of 380-580mm (15"-23") adjustable between 380-480mm (15"- 19") and 480-580mm (19"-23").
- Backrest heights of 560-660mm in 25mm increments
- Cantilever Flip-Back Arm, width, height and angle adjustable with Desk armpad with Startex covered visco foam
- Stealth all positioning headrest
- Ultra Low Maxx CG power tilt (0 to 50°) seat heights of 420-480mm
- Powered Recline 90° to 170° with Extended Shear Reduction (ESR)
- Angle adjustable legrest hanger (swinging away and detachable) Vari F with Fixed footplates
- DX2 control system with REM 550 with LCD screen
- Centre (Driving Wheels) 3.00 by 8.00 Grey Puncture Proof (Solid Inlay)
- Selection of Vibrant Chassis Shroud Colours
- 10km/h (4mph) 4 Pole Motors and 73Ah Batteries
- G-Trac
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



**Standard Base Price : £9,950**

x Please make your selection with a cross in the box ◆ Standard feature ○ Free of Charge

Code	Description	DEALER
ARO0080	TDX SP2 ULTRA	£9,950
ARZ0080	TDX SP2 ULTRA Special Basic Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM	POA
<b>SEATING TYPE : ULTRA LOW MAXX</b>		
ARO2020	X ULTRA LOW MAXX Seating Unit	◆

SEAT DIMENSIONS												
		Seat Frame Width Range										
		ARO2021 405MM-505MM (INCREMENTS OF 25MM)/ 16"-20"					ARO2022 505MM-610MM (INCREMENTS OF 25MM)/ 20"-24"					
		Seat width setting for adjustable seat frame										
		ARO2025 405MM/ 16"	ARO2026 430MM/ 17"	ARO2027 455MM / 18"	ARO2028 480MM / 19"	ARO2029 505MM / 20"	ARO2029 505MM / 20"	ARO2030 530MM / 21"	ARO2031 555MM / 22"	ARO2032 580MM / 23"	ARO2033 610MM / 24"	
Seat Frame Depth Range	19"	ARO2034 380 MM / 15"	○	○	○	○	○	○				
		ARO2035 405MM / 16"	○	○	○	○	○	○				
		ARO2036 430 MM / 17"	○	○	○	○	○	○				
	23"	ARO2037 455 MM / 18"	○	○	○	○	○	○	○	○	○	○
		ARO2038 480MM / 19"	○	○	○	○	○	○	○	○	○	○
		ARO2038 480MM / 19"	○	○	○	○	○	○	○	○	○	○
		ARO2039 505 MM / 20"	○	○	○	○	○	○	○	○	○	○
		ARO2040 530MM / 21"	○	○	○	○	○	○	○	○	○	○
		ARO2041 555 MM / 22"			○	○	○	○	○	○	○	○
		ARO2042 580 MM / 23"			○	○	○	○	○	○	○	○

NOTE: Seat Depths 380 mm (15") and 405mm (16") are not possible with Centre mounted legrest

**SEAT CUSHION (sizes according to seat frame settings)**

ARO2144	Matrix Vi Cushion (Mesh fabric)		£335
ARO2146	Matrix Flo-tech Contour Visco	Not Available in Seat depth of 580mm (2042) or 555mm (2041) in combination with seat widths of 405mm (2025), 430mm (2026), 455mm (2027) 480mm (2028)	£92
ARO2147	Matrix Libra Cushion, 90mm thick	Not available with Seat Width 530mm with 410mm seat Depth	£319
ARO2140	Solid seat Base, For alternative cushion options (Vicair, Libra etc) please see refer to our seating/positioning product ranges.		◆




**BACKREST SUPPORT OPTIONS (material height, see below for overall back height)**

ARO2050	High Back, height 410 mm (Mesh fabric)	Not possible with backrest height of 635mm and 660mm	◆
ARO2051	High Back, height 460 mm (Mesh fabric)	Not available with backrest height of 560mm and 585mm	○
ARO2174	Matrix Posture Back, height 400 mm (Mesh fabric max dimension 455 mm wide)	Not available with seat widths of 480mm (19")- 610mm (24"), backrest height of 610mm, 635mm and 660mm	£406
ARO2175	Matrix Posture Back, height 500 mm (Mesh fabric max dimension 455 mm wide)	Not available with seat widths of 480mm (19")- 610mm (24"), backrest height of 560mm, 585mm, 610mm	£406
ARO2177	Matrix Elite TR Back, height 500 mm (Mesh fabric)	Not available with backrest height of 560mm, 585mm, 610mm, 635mm	£406
ARO2178	Matrix Elite Deep Back, height 500 mm (Mesh fabric)	Not available with backrest height of 560mm, 585mm, 610mm, 635mm	£450





**BACKREST HEIGHT: Measured from seat pan to top of the back shell**

**NOTE: important for the adjustment of the ESR recline**




ARO2070	560 mm (22")	◆
ARO2071	585 mm (23")	○
ARO2072	610 mm (24")	○
ARO2073	635 mm (25")	○
ARO2074	660 mm (26")	○

<b>ARMREST</b>			
ARO2095		Cantilever Flip-Back Arm, width, height and angle adjustable (240 - 330 mm)	◆
ARO2096		Cantilever Flip-Back Arm, width, height and angle adjustable (320 - 405 mm)	○
ARO2097		Reclining Arm, Flip Back, adjustable height (240 - 330 mm)	○
ARO2098		Reclining Arm, Flip Back, adjustable height (330 - 405 mm)	○
ARO0690		Without Armrest	○
<b>ARMPAD</b>			
<b>Left Hand</b>			
ARO2012		Desk arm pad with startex cover visco foam (255mm)	◆
ARO2013		Full length arm pad with startex cover visco foam (355mm)	○
ARO2014		Ergonomic arm trough	£80
ARO2015		Without LH arm pad	○
<b>Right Hand</b>			
ARO2016		Desk arm pad with startex cover visco foam (255mm)	◆
ARO2017		Full length arm pad with startex cover visco foam (355mm)	○
ARO2018		Ergonomic arm trough	£80
ARO2019		Without RH arm pad	○
<b>ARM PAD ACCESSORIES &amp; MOUNTING OPTIONS</b> (in combination with arm trough only); to be mounted on the opposite side of remote!			
<b>Left Hand</b>			
ARO2080		Flat Hand Pad Left	£80
ARO2081		Ergonomic Mounting Hardware to mount Otto Bock (EMO) hand pad options Left	£35
<b>Right Hand</b>			
ARO2083		Flat Hand Pad Right	£80
ARO2084		Ergonomic Mounting Hardware to mount Otto Bock (EMO) hand pad options Right	£35
<b>MULTI-AXIS UPPER EXTREMITY SUPPORTS</b>			
<b>Left Hand</b>			
ARO2082		Multi Axis Upper Extremity Support (MACES) Left,raises arm pad by 43 mm	£200
<b>Right Hand</b>			
ARO2085		Multi Axis Upper Extremity Support (MACES) Right,raises arm pad by 43 mm	£200
<b>Elbow Support LH</b>			
ARO2316		Swing away Elbow Block LH, (pad size 100 x 100 mm)	Not available with ergonomic arm trough LH £375
ARO2317		Swing away Elbow Block LH, (pad size 75 x 250 mm)	Not available with ergonomic arm trough LH £375
<b>Elbow Support RH</b>			
ARO2318		Swing away Elbow Block RH, (pad size 100 x 100 mm)	Not available with ergonomic arm trough RH £375
ARO2319		Swing away Elbow Block RH, (pad size 75 x 250 mm)	Not available with ergonomic arm trough RH £375
<b>HIP SUPPORTS LH (for flip up and following armrests only)</b>			
ARO2184		Multi adjustable hip support LH, pad size SMALL 140 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
ARO2186		Multi adjustable hip support LH, pad size MEDIUM 230 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
ARO2188		Multi adjustable hip support LH, pad size LARGE 330 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
<b>HIP SUPPORTS RH (for flip up and following armrests only)</b>			
ARO2185		Multi adjustable hip support RH, pad size SMALL 140 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
ARO2187		Multi adjustable hip support RH, pad size MEDIUM 230 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
ARO2189		Multi adjustable hip support RH, pad size LARGE 330 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
<b>SEAT ACCESSORIES</b>			
<b>HEADREST AND MOUNTING HARDWARE</b>			
ARO2092		Onyx Auto Style headrest	£90
ARO2093		Comfort Plus with 255 mm Pad	£90
ARO2094		All Positioning Headrest Large	◆
 <b>Onyx Auto Style headrest</b>  <b>Comfort Plus with 255mm pad</b>  <b>All Positioning Headrest Large</b>			

<b>LATERAL SUPPORTS LH</b>				
ARO2043		Lateral Support, Quick Release, pad size (145 mm wide x 95 mm high) LH	£245	
ARO2045		Lateral Support, Quick Release, pad size (165 mm wide x 110 mm high) LH	£290	
<b>LATERAL SUPPORTS RH</b>				
ARO2044		Lateral Support, Quick Release, pad size (145 mm wide x 95 mm high) RH	£245	
ARO2046		Lateral Support, Quick Release, pad size (165 mm wide x 110 mm high) RH	£290	
<b>LATERAL SUPPORT ACCESSORIES</b>				
ARO2049		Additional telescopic link for lateral support, 25 mm, Pair	£115	
<b>SEAT ADJUSTMENTS</b>				
ARO2052		Centre of gravity Tilt (0 to 50°)	◆	
ARO2053		300 mm seat lifter (increased STF Height by 4 mm) with centre of gravity tilt (0 to 45°)	Not Possible with seat-to -floor height of 480mm £1,750	
<b>SEAT HEIGHT SETTING</b>				
ARO0337		Seat-to-Floor Height 440 mm	○	
ARO0331		Seat-to-Floor Height 460 mm	○	
ARO0339		Seat-to-Floor Height 480 mm	○	
<b>ANTERIOR and POSTERIOR TILT WEDGE CAN BE EASILY RETRO FITTED OR REMOVED (additional wedge situated between seat frame a Tilt/Lifter unit)</b>				
ARO2057		Tilt Standard Mount (no anterior or posterior tilt)	◆	
ARO2054		5° Fixed Anterior Assist Mount (will reduce max tilt angle by 5°) <b>Front STF height as defined above</b>	£295	
ARO2055		10° Fixed Anterior Assist Mount (will reduce max tilt angle by 10°) <b>Front STF height as defined above</b>	£295	
ARO2056		5° Pre-Tilt (will increase max tilt angle by 5°) <b>Front STF height increased by 60 mm over chosen seat height setting</b>	£115	
<b>ACTUAL SEAT TO GROUND HEIGHTS WITH ANTERIOR AND POSTERIOR OPTIONS</b>				
	<b>Seat To Ground Height</b>	<b>5° Anterior Tilt ARO2054</b>	<b>10° Anterior Tilt ARO2055</b>	<b>5° Pre Tilt ARO2056</b>
	440 mm	440 mm	440 mm	500 mm
	460 mm	460 mm	460 mm	520mm
	480 mm	480mm	480 mm	540 mm
<b>BACKREST RECLINER</b>				
ARO2060		Powered Recline 90 to 170° with Extended Shear Reduction (ESR), Standard Mount	◆	
ARO2063		Powered Recline, 8° Precline mount, 82° to 162° with Extended Shear Reduction (ESR)	£25	
ARO2066		Powered Recline, 30° Precline mount, 60° to 140° with Extended Shear Reduction (ESR)	£25	
<b>PLEASE CHOOSE EITHER CENTRE MOUNTED OR SIDE MOUNTED LEGRESTS</b>				
<b>CENTRE MOUNTED LEGRESTS</b>				
<b>FOOTREST HANGER TYPE (INCLUDING CALF PAD AND ANGLE ADJ FOOTPLATE)</b>				
ARO0586		Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 390 - 510mm	£233	
ARO0589		Hanger, centre mounted, length compensating, powered, LNX (0° - 90°): knee to heel length 390 - 510mm	£1,300	
<b>FOOTPLATES</b>				
ARO0559		Angle & width adjustable footplates (240 mm depth)	○	
ARO0516		Angle & width adjustable footplates (290 mm depth)	○	
ARO0517		Angle adjustable foot platform (240 mm depth, 305 mm width)	£45	
<b>SIDE MOUNTED LEGRESTS</b>				
<b>INDIVIDUAL FOOTREST HANGER TYPE; LH</b>				
ARO0536		<b>VARI F</b> : Angle adjustable legrest hanger (swing away and detachable) (adjustable via tools 0° - 70°) : knee to heel length 290 - 460mm	◆	
ARO0537		<b>VARI A</b> : Manually angle adjustable legrest hanger (swing away and detachable) (adjustable via hand wheel / knob 0° - 70°) - complete with calf pad support : knee to heel length 290 - 460mm	Not possible with Fixed footplate LH £89	

ARO0538		<b>ADM Manual</b> : Manually angle adjustable (0° - 80°) legrest hanger (swing away and detachable) (adjustable via quick release lever) - complete with calf pad support : knee to heel length 290 - 460mm	Not possible with Fixed footplate LH	<b>£158</b>
ARO0539		<b>ADE Powered</b> : Electrically angle adjustable (0° - 80°) legrest hanger (swing away and detachable) (powered via joystick) - complete with calf pad support : knee to heel length 290 - 460mm	Not possible with Fixed footplate LH	<b>£478</b>
ARO0572		Pivot Plus Manual Elevating Legrests LH	Not possible with Fixed and Angle Adjustable footplate LH	<b>£520</b>
ARO0577		Pivot Plus Power Elevating Legrests LH	Not possible with Fixed and Angle Adjustable footplate LH	<b>£850</b>
ARO0535		Without Hanger LH and footplates		○
<b>INDIVIDUAL FOOTREST HANGER TYPE; RH</b>				
ARO0546		<b>VARI F</b> : Angle adjustable legrest hanger (swing away and detachable) (adjustable via tools 0° - 70°) : knee to heel length 290 - 460mm		◆
ARO0547		<b>VARI A</b> : Manually angle adjustable legrest hanger (swing away and detachable) (adjustable via hand wheel / knob 0° - 70°) - complete with calf pad support : knee to heel length 290 - 460mm	Not Possible with Fixed Footplate	<b>£89</b>
ARO0548		<b>ADM Manual</b> : Manually angle adjustable (0° - 80°) legrest hanger (swing away and detachable) (adjustable via quick release lever) - complete with calf pad support : knee to heel length 290 - 460mm	Not Possible with Fixed Footplate	<b>£158</b>
ARO0549		<b>ADE Powered</b> : Electrically angle adjustable (0° - 80°) legrest hanger (swing away and detachable) (powered via joystick) - complete with calf pad support : knee to heel length 290 - 460mm	Not Possible with Fixed Footplate	<b>£478</b>
ARO0579		Pivot Plus Manual Elevating Legrests RH	Not Possible with Fixed Footplate or Angle adjustable footplate	<b>£520</b>
ARO0584		Pivot Plus Power Elevating Legrests RH	Not Possible with Fixed Footplate or Angle adjustable footplate	<b>£850</b>
ARO0545		Without Hanger LH and footplate		○
<b>INDIVIDUAL FOOTPLATE TYPE; LH</b>				
ARO0505		Fixed Footplate LH		◆
ARO0509		Angle Adjustable Footplate LH		<b>£52</b>
ARO0511		Angle and depth adjustable Footplate LH		○
<b>INDIVIDUAL FOOTPLATE TYPE; RH</b>				
ARO0502		Fixed Footplate RH		◆
ARO0506		Angle Adjustable Footplate RH		<b>£52</b>
ARO0508		Angle and depth adjustable Footplate RH		○
<b>ACS2 CONTROLS</b>				
<b>USER CONTROL + BRACKET</b>				
ARO1617		REM 550 with LCD screen	Not possible with Chin remote	◆
ARO1601		ACS REM24SD	Not possible with Quad link swing away	○
ARO1614		G90A	Not possible with Chin remote or Quad link swing away	<b>£147</b>
ARO1615		G91S Includes bracketry, select option 1207 or 1208	Not available with swing away bracketry High end or Quad link swing away	<b>£866</b>
		Head-Remote, ultrasonic proximity sensor(only available as special) Please contact UK Specials Dept		<b>POA</b>
 <b>REM 550</b>  <b>REM 24SD</b>  <b>G90 A</b>  <b>G91 S</b>				

ARO1207		Standard Bracket RH		○
ARO1208		Standard Bracket LH		○
ARO1707		Swing Away + height adjust.RH High End		£105
ARO1708		Swing Away + height adjust. LH High End		£105
ARO2075		Swing Away + height / width adjust. Mount on RH, Quad link design		£165
ARO2076		Swing Away + height / width adjust. Mount on LH, Quad link design		£165



**ATTENDANT REMOTE + BRACKET Note: Not possible in combination with IDC**

ARO1618		REM 550 with LCD screen	Not possible with compact dual control or chin remote.	£212
ARO1602		ACS REM24SD	Not possible with compact dual control or chin remote.	£158
ARO1616		G90A	Not possible with Chin remote or compact dual control	£158
ARO1207		Standard Bracket RH		○
ARO1208		Standard Bracket LH		○

**TRAY & CENTRE CONTROLS**

ARO1608		A C S Centre Control (REM 41)	Not possible with Chin remote	£168
ARO1607		A C S Tray Control c/w Acrylic Tray & Actuator Functionality - DX REM 41	Not possible with Chin remote	£268
ARO1289		Swing Away RH		○
ARO1288		Swing Away LH		○


**CHIN CONTROL SIMPLE**

ARO1719		Chin Remote with standard Joystick	Only available with Attendant Rem24SD as special	£1,646
ARO1289		New Manual Swing Away Bracket RH	Not possible with RH IDC	○
ARO1288		New Manual Swing Away Bracket LH	Not possible with LH IDC	○
ARO1287		New Swing Away Bracket RH, powered	Not possible with RH IDC	£1,095
ARO1286		New Swing Away Bracket LH, powered	Not possible with LH IDC	£1,095


**DUAL CONTROL + BRACKET**

ARO1720		Compact Dual Control		£377
ARO1207		Standard Bracket RH		○
ARO1208		Standard Bracket LH		○

**Dual control**



**IDC Dual control**




**INVACARE IDC DUAL CONTROL + BRACKET (Only available with ACS)**

ARO1645		Invacare IDC Dual control		£542
ARO1221		Detachable Bracket RH		○
ARO1222		Detachable Bracket LH		○

**DRIVE UNIT OPTIONS**

**FRAME COLOURS (ALL COLOURED PARTS ARE DETACHABLE AND CHANGEABLE PLASTIC SHROUDS)**

ARO0228		Wet Black		○
ARO0212		Electric Red		○
ARO0257		Jade Green		○
ARO0236		Electric Blue		○
ARO0237		Electric Silver		○
ARO0238		Arctic White		○
ARO0205		Transparent Purple		○
SWO0001		Carbon Effect (Extended lead time)		○
SWO0002		True Flames Blue (Extended lead time)		○
SWO0003		True Flames Green (Extended lead time)		○

CENTRE WHEELS (DRIVING WHEELS)			
ARO0750		14" diameter x 3" Grey Puncture Proof (Solid Inlay) Drive wheels	◆
ARO0773		14" diameter x 3" Grey Puncture Protected (Sealed) Drive wheels	○
ARO0771		14" diameter x 3" Grey Pneumatic	○
PERFORMANCE OPTIONS			
MOTORTYPE			
ARO1282		6 km/h (4mph)	○
ARO1285		10 km/h (6mph) (includes light system and indicators - please select option ARO1437)	◆
BATTERY			
ARO1260	X	73Ah	◆
G-TRACK			
ARO1215	X	 <b>G-trac</b> ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain.	◆
CHARGER			
ARO1343	X	Charger 8A (200 - 250V)	◆
SAFETY OPTIONS			
LIGHT			
ARO1437		Light system with Indicators & LED Headlights	£286
ARO1254		Without Light system	◆
SAFETY ACCESSORIES			
TRANSPORTATION KIT			
ARO1450		<b>Preparation</b> for DAHL Docking station: Docking station for car implementation, please order docking station separately (Adequate headrest is compulsory, see options above)	£430
ARO1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
HEADRESTS ADEQUATE FOR PERSON TRANSPORTATION			
		INVACARE Headrest (Please select a headrest from this form)	◆
FIXATION BELT			
ARO2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	◆
REARVIEW MIRROR (mounted opposite side to controller)			
ARO1415		RH	○
ARO1339		LH	○
ACCESSORIES			
TRAY - mounted opposite side to controller			
ARO1404		Tray Swing Away RH	Not available with 610mm (24") seat width £143
ARO1405		Tray Swing Away LH	Not available with 610mm (24") seat width £143
ARO2310		Flip and fold away 1/2 lap tray LH	£605
ARO2311		Flip and fold away 1/2 lap tray RH	£605
KLiCKfix BAG AND BRACKETS			
ARO1480		HandyZip bag, removable	£51
ARO1478		Quick release bracket RH	Not possible with Rh swing away control ○
ARO1479		Quick release bracket LH	Not possible with Lh swing away controls ○
<b>TOTAL PRICE</b>			



**Special Instructions - Please let us know your preferred contact method:**

**Name:**

**Email:**

**Fax:**

**Description of special request (please include relevant option code or part number):**