

SHOPPRIDER



Whisper

Model No. **S787NA**

Special Features

- Easy to dismantle
- Two colour options
- Heaviest part 17kgs
- Flip up arms
- Adjustable length



Specifications



Model No.*	S787NA	Battery Size 12v x 2	12 Ah
Class	2	Range**	Up to 20 km / 12 miles
Overall Length	940-1020 mm / 37"- 40"	Overall Weight	39 kg / 86 lbs
Overall Width	430 mm / 17"	Maximum Speed	6 kmph / 4mph
Overall Height	870 mm / 34"	Wheel / Tyre Size	175 & 200 mm / 7 & 8"
Seat Width	400 - 500 mm / 16"- 20"	Turning Radius	1000 - 1200 mm / 42"- 50"
Seat Height (to ground)	460 - 500 mm / 18"- 20"	Battery Charger	1.2 amp off board
Armrest Height	150 - 200 mm / 6"- 8"	Ground Clearance***	60 mm / 2.5"
Backrest Height from Seat (inc headrest)	300 - 320 mm / 12"- 13"	Colours	Blue
Seat to Footrest/Floorpan Height	330 - 380 mm / 13"- 15"		
Safe Working Slope	6° / 10%		
Maximum User Weight	115 kg / 18 st		

* Lower limb amputees should refer to user manual before use. ** Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel. *** Measurement to underside of chassis