

# Torino

Model No. **S778XLSBN**

**SHOPPRIDER**



## Special Features

- Detachable fold-down lockable 4 position swivel Captain seat with forward, back and height adjustments
- Flip-up arms with width, and angle adjustment
- Full suspension
- Black non-marking tyred wheels
- Free wheel release handle
- Back indicator light cover
- Delta handle
- Five position tiller handle

### Optional Extras

- Rear Basket
- Crutch and walking stick holder
- Canopy

## Specifications



Model No.*	<b>S778XLSBN</b>	Maximum User Weight	<b>220 kg / 35 st</b>
Class	<b>3/C</b>	Battery Size 12v x 2	<b>75 Ah</b>
Overall Length	<b>1360 mm / 53"</b>	Range**	<b>Up to 55 km / 35 miles</b>
Overall Width	<b>680 mm / 27.5"</b>	Overall Weight	<b>140 kgs / 308 lbs</b>
Overall Height	<b>1270 mm / 50"</b>	Maximum Speed	<b>6 &amp; 12 kmph / 4 &amp; 8 mph</b>
Seat Width	<b>500 mm / 20"</b>	Wheel / Tyre Size (Front)	<b>300 mm / 12"</b>
Seat Depth	<b>440 mm / 17"</b>	(Back)	<b>360 mm / 14"</b>
Seat Height (to ground)	<b>630 mm-705 mm / 25"-28"</b>	Turning Radius	<b>1250 mm / 49"</b>
Armrest Height	<b>210 mm / 8"</b>	Battery Charger	<b>8 amp off board</b>
Backrest Height from Seat (inc headrest)	<b>660-790 mm / 26"- 31"</b>	Ground Clearance	<b>150 mm / 6"</b>
Seat to Footrest/Floorpan Height	<b>430-505 mm / 17"-20"</b>		
Safe Working Slope	<b>8° / 14%</b>		

\* Lower limb amputees should refer to user manual before use. \*\* Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.