



A Light, Portable
Mid-range Scooter



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

Pride®
Mobility Products Ltd.
— LIVE YOUR BEST® —



Tel 0333 305 8274

www.pride-mobilityscooters.co.uk

admin@pride-mobilityscooters.co.uk

Features and Specification

Introduction...

Delivering reliable performance and featuring a practical yet stylish design. The Apex Sprint is extremely portable, durable and comfortable. The Sprint fits in almost any car boot and offers plenty of legroom for all users. In addition the Apex Sprint offers a stable secure base and pneumatic tyres, making it an excellent all round scooter.

Quick and easily detachable

The Apex Sprint's unique chassis enables quick and easy detachment. Allowing it to be disassembled into 5 pieces, which can be safely stored in most small spaces.



Unique Delta Tiller Design

Featuring a wrap-around tiller design. The Apex Sprint enables fine adjustments that best suit your needs.



True Blue



Champagne Orange

Available Colours

The Apex Sprint comes in two optional colours:

*images used for illustration purposes only
(Version 1: April 2015)

Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
21 Stone/6lbs. (136 kg)	4 Mph (6.4 km/h)	6 cm (2.36 inch's)	127 cm (50 inch's)	114 cm (44.9 inch's)	55 cm (21.7inch's)	Front: 19cm Rear: 23 cm	8.4 Stone 53.5 k.g.	3.4 Stone (Front Sec)	Folded Seat 12.6 k.g	x2 35 AH deep cycle batt

Warranty

Two-year limited warranty on frame
One-and-a-half-year limited warranty on Controllers
One-year limited warranty Electronic Components, battery and Charger.

Notes

Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer.

Tel 0333 305 8274

www.pride-mobilityscooters.co.uk

admin@pride-mobilityscooters.co.uk

Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST —