

GO-GO[®]
ELITE TRAVELLER PLUS[®]

THE QUICKEST WAY
TO TRAVEL.

GO-GO ELITE TRAVELLER[®] PLUS

The Go-Go Elite Traveller[®] Plus was designed to bring travel scooters to a greater range of people than ever before. An increase in length and width, and a weight capacity of 147 kg (23 Stone 3 lbs.) combine with wraparound delta tiller to allow both the larger individual and those with limited dexterity to experience the benefits of a specific travel mobility product for the first time. Feather-touch disassembly makes the Go-Go Elite Traveller Plus the easiest travel scooter to take along anywhere.



(Images are for illustration purposes only)



Tel. 0333 305 8274

www.pride-mobilityscooters.co.uk

admin@pride-mobilityscooters.co.uk

Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST[®] —

Go-Go Elite Traveller® Plus	3-wheel	4-wheel
Weight capacity:	147 kg (23 Stone 3 lbs.)	147 kg (23 Stone 3 lbs.)
Maximum speed ¹ :	Up to 6.4 km/h (4 mph)	Up to 6.4 km/h (4 mph)
Ground clearance ² :	7.62 cm (3") at midframe	7.62 cm (3") at midframe
Turning radius ² :	87.6 cm (34.5")	116.2 cm (45.75")
Overall length ² :	101.6 cm (40")	102.25 cm (40.25")
Overall width ² :	53.7 cm (21.2")	53.7 cm (21.2")
Tyre size:		
Front:	5 cm x 20.32 cm (2" x 8") solid	5 cm x 20.32 cm (2" x 8") solid
Rear:	7.62 cm x 22.86 cm (3" x 9") solid	7.62 cm x 22.86 cm (3" x 9") solid
Maximum allowable incline:	10.5% (6°)	10.5% (6°)
Maximum angle capacity:	17.6% (10°)	17.6% (10°)
Maximum obstacle climbing ability:	6.5 cm (2.5")	6.5 cm (2.5")
Total weight w/o battery pack ⁴ :	39 kg (86 lbs.)	39.2 kg (86.5 lbs.)
Heaviest piece when disassembled:	14.7 kg (32.5 lbs.) (front section)	16 kg (35 lbs.) (front section)
Drive system:	24-volt DC motor, sealed mini-transaxle	
Dual braking system:	Electronic, regenerative, and electromechanical	
Standard seating:	Type: Foldable Weight: 10.2 kg (22.5 lbs.) Material: Black vinyl Dimensions: 43.18 cm (17") overall width 43.18 cm (17") effective depth	
Battery requirements ³ :	Two 12-volt, deep cycle (standard) Size: 2 x 17 AH (standard) Weight ⁵ : 17 AH: 13.4 kg (29.5 lbs.) battery pack	
Battery charger:	Off-board, 2-amp	

Standard colours:



Blue



Red

2 sets of exciting colours with every scooter!

Features

- Feathertouch disassembly with one hand
- Exclusive black, lowprofile, non scuffing tyres with new trispoke mags
- Delta tiller with wraparound handles
- 23 cm (9") rear wheels
- 147 kg (23 Stone 3 lbs.) weight capacity
- Autoconnecting front to rear cable
- Includes 2 sets of easily changeable coloured panels Red and Blue
- Front framemounted seat post offers maximum stability
- All new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Microprocessorbased controller offers optimal power management and added safety features
- Front tillermounted basket
- Removable footboard basket for additional storage (4wheel version only)
- Convenient offboard dual voltage charger can charge battery pack onboard or offboard

Options

- Cup holder
- Rear basket
- Crutch holder



Feathertouch disassembly

¹ Varies with user weight, terrain type, battery ampour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of + 10%, 5%.

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).

³ AGM or gel cell type required.

⁴ Includes standard seat.

⁵ Battery weight may vary based on manufacturer.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 series EN12184 standards. All specifications subject to change without notice.