

4-WHEEL SCOOTER



THE ULTIMATE MID-SIZE LUXURY MOBILITY SCOOTER

The **Apex Finesse** Mobility Scooter features Comfort-Trac front and rear independent suspension for high performance. Offering a sleek, modern design and a full complement of standard features, including a delta tiller with wraparound handles and a luggage compartment light, the Apex Finesse is an excellent mobility scooter option.



Front Suspension



Rear Suspension



Delta tiller with wraparound handles

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

www.pride-mobilityscooters.co.uk

SPECIFICATIONS


APEX FINESSE

Weight capacity:	28 Stone 8 lbs. (181.4 kg)
Maximum speed ¹ :	Up to 4 mph (6.4 km/h)
Ground clearance ² :	5.6 cm (2.2") (at motor)
Turning radius ² :	152.1 cm (59.9")
Overall length ² :	119.4 cm (47")
Overall width ² :	56.5 cm (22.2")
Seattofloor height:	45.1 cm (17.7")
Front tyres ² :	25.4 cm (10") (pneumatic)
Rear tyres ² :	25.4 cm (10") (pneumatic)
Maximum allowable incline:	6°
Total weight without batteries ⁴ :	70.9 kg (156.2 lbs.)
Heaviest piece when disassembled:	35.2 kg (77.6 lbs.) (front section)
Standard seating:	Type: Folded Moldable Plastic Weight: 12.6 kg (27.8 lbs.) Material: Black Vinyl Dimensions: 43.18 cm x 43.18 cm (17" x 17")
Motor output:	24-volt, DC motor
Brakes:	Regenerative and electromechanical
Battery requirements ⁵ :	(2) 12-volt, deep cycle Size: (2) 35 AH Weight ³ : (35 AH) 9.7 kg (21.4 lbs.) each
Battery charger:	5-amp, off-board

Colour:



Blue



Red

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).

³ Battery weight may vary based on manufacturer.

⁴ Includes standard seat weight.

⁵ AGM or gel cell type required.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

DEALER INFO HERE

UK_APEX FINESSE_SS_3-25-15

FEATURES

- CTS Suspension (Comfort-Trac front and rear independent suspension)
- Sleek, modern design
- Exclusive black, nonscuffing tyres
- Modular design for easy serviceability
- Delta tiller with wraparound handles
- Easily accessible charger port in tiller
- Full LED lighting package that includes kerb and rear lights and directional signals.
- Courtesy light under tiller
- Storage compartment below tiller
- Enhanced flushmounted horn buttons
- Auto-connecting front-to-rear cable

OPTIONS

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Walker holder
- Saddlebag
- Rearview mirror



Rear LED lights

Pride
Mobility Products BV
— LIVE YOUR BEST —

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

Tel. 0333 305 8274

www.pride-mobilityscooters.co.uk

admin@pride-mobilityscooters.co.uk



British
Healthcare
Trades
Association

