

<b>Invoice to:</b>	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Order date:	
Purchase order:	

<b>Deliver to:</b>	
Name:	
Street:	
Town:	
Country:	
Tel:	
Marked for:	

<b>Standard Features</b>	<ul style="list-style-type: none"> <li>● Angle adj. footplate or central legrest</li> <li>● VR2 90Ah controller</li> <li>● Parallel swing away mount</li> <li>● 6KPH (4-pole) motors</li> </ul>	<ul style="list-style-type: none"> <li>● Solid drive wheels</li> <li>● Solid castors</li> <li>● 60Ah batteries</li> <li>● Std curved headrest</li> </ul>
--------------------------	---	--

Fax: +44 (0)845 605 66 89

E-mail: orders@sunmed.co.uk

Spare Parts: <http://parts.sunrisemedical.eu>

**JIVE R2: £5351**

**JIVE F: £6175**

**JIVE M: £6644**

**Max. user weight: 160kg (25 stone)**

**SEAT DIMENSION JAY COMFORT**

Seat frame width is 41 cm. The Seat width can be adjusted between 40 and 56 cm via the armrests, based on the seat/back cushion configuration. The Seat depth can be adjusted between 41 and 55 cm via backrest movement  
 Seat dimensions also reflect aluminium seat plate option. Seat plate option to be selected at the Jay Comfort Cushion section (JIV020109).

**Cushion width**

		inches		15.5	16.5	17	18	19	20	20.5	21	22
		cm		40	42	44	46	48	50	52	54	56
<b>Seat depth setting</b>	48 cm depth seat cushion	16	41	○	○	○	○	○	○	○	○	○
				JIV020001	JIV020021	JIV020005	JIV020025	JIV020009	JIV020029	JIV020013	JIV020033	JIV020017
	55 cm depth seat cushion	18	46	○	○	○	○	○	○	○	○	○
				JIV020002	JIV020022	JIV020006	JIV020026	JIV020010	JIV020030	JIV020014	JIV020034	JIV020018
	20	20	51	○	○	○	○	○	○	○	○	○
				JIV020003	JIV020023	JIV020007	JIV020027	JIV020011	JIV020031	JIV020015	JIV020035	JIV020019
	21.5	55	○	○	○	○	○	○	○	○	○	
				JIV020004	JIV020024	JIV020008	JIV020028	JIV020012	JIV020032	JIV020016	JIV020036	JIV020020

**SEAT HEIGHT & DRIVE WHEEL POSITION**

**Seat height**



inches	16.5	17	17.5	18	18.5	19
cm	42	43	44.5	46	47	48
<b>RWD</b>		○		○(3)		
JIV010011		JIVH43		JIVH46		
<b>FWD</b>		○(1)	○	○	○	○
JIV010002		JIVH43	JIVH44,5	JIVH46	JIVH47	JIVH48
<b>MWD</b>	○ (2)	○	○	○	○	○
JIV010003	JIVH42	JIVH43	JIVH44,5	JIVH46	JIVH47	JIVH48

Seat height measured floor to seat pan at 0° tilt

**When powered lift/tilt options are selected achievable seat height is as follows:**

- FWD: Tilt, only: 44.5cm, Lift & Tilt: 46cm
- MWD: Tilt: 43cm, Lift & Tilt: 44.5cm
- RWD: Tilt only: 43cm, Lift & Tilt: 44.5cm

JIV090022 ○ Seat height extension - Adds 3cm to seat height selected (3) 30mm seat height kit to be selected (JIV090022)

(1) not available with power modules or 0° seat angle setting with no module  
 (2) not available with filler module in 0° and 3° seat angle setting (JIV030001)

**SEAT TO BASE INTERFACE**

JIV030001	<input type="radio"/>	Fixed tilt (filler module) <span style="float: right;">(tool adjustable seat angle)</span>	
	<input type="radio"/>	0° - (Jive F / Jive M only)	<input type="radio"/>
	<input type="radio"/>	3°	<input type="radio"/>
	<input type="radio"/>	6°	<input type="radio"/>
	<input type="radio"/>	9°	
<b>Power seating modules</b>			
JIV030002	<input type="checkbox"/>	30° Power COG tilt	592
JIV030003	<input type="checkbox"/>	50° Power COG tilt <span style="float: right;">(Jive MWD Only)</span>	631
JIV030004	<input type="checkbox"/>	-5° to 25° Powered COG tilt <span style="float: right;">(3cm added seat to floor height at 0°)</span>	695
JIV030006	<input type="checkbox"/>	Power seat lift and tilt combined <span style="float: right;">(30° tilt and 30cm lift)</span>	1589

**JAY CUSHION SIZE**

**SEAT WIDTH**

		Width Coding	JL40	JL42	JL44	JL46	JL48	JL50	JL56
Depth Coding	cm	40	42	44	46	48	50	56	
JL__40	40	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>			
JL__42	42	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>			
JL__44	44	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>			
JL__46	46	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
JL__50	50				<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
JL__56	56							<input type="radio"/>	<input type="radio"/>

**Size Configuration - COMPLETE PRODUCT**

**Conversion Table**

cm	40.0	42.0	44.0	46.0	48.0	50.0	56.0
inches	15.6	16.4	17.2	18.0	18.8	19.5	21.9

(1) Not available with Soft Combi P or Easy Fluid Cushion

(2) Not available with JAY Lite standard version

**JAY CUSHION OPTIONS**

JAY020004	<input type="checkbox"/>	Jay Lite Standard version <span style="float: right;">(Micro climat cover, Aluminium seat plate to be selected)</span>	195
JAY020005	<input type="checkbox"/>	Jay Easy Fluid with Flat Base <span style="float: right;">(Incontinence cover, Aluminium seat plate to be selected)</span>	280
JAY020001	<input type="checkbox"/>	Jay Easy Visco with Flat Base <span style="float: right;">(Incontinence cover, Aluminium seat plate to be selected)</span>	120
JAY020003	<input type="checkbox"/>	Jay Soft Combi P <span style="float: right;">(Incontinence cover, Aluminium seat plate to be selected)</span>	70

**JAY COMFORT SEAT CUSHION OPTIONS**

JIV020109	<input type="checkbox"/>	Aluminium seat pan - no cushion required <span style="float: right;">(Includes omit cushion option JIV020209)</span>	-106
JIV020101	<input type="radio"/>	Comfort seat pan, without threshold <span style="float: right;">(thickness 2cm)</span>	
JIV020102	<input type="radio"/>	Comfort seat pan, with 5cm (2") threshold <span style="float: right;">(thickness 2cm)</span>	
JIV020201	<input type="radio"/>	Standard cushion - stretched rubber upholstery <span style="float: right;">(thickness 8cm)</span>	
JIV020202	<input type="radio"/>	Optimum comfort - stretched rubber upholstery <span style="float: right;">(anatomic formed) (thickness 9cm)</span>	



JIV020101



JIV020102



JIV020201



JIV020202

**JAY COMFORT BACK CUSHION OPTION**

(57cm back height - Adj. up to 70cm)

JIV040001	<input type="radio"/>	Shallow contour (standard velour) <span style="float: right;">(hard back)</span>	
JIV040002	<input type="checkbox"/>	Shallow contour (standard velour) <span style="float: right;">(vertical elastic straps)</span>	101
JIV040003	<input type="checkbox"/>	Shallow contour (standard velour) <span style="float: right;">(horizontal tension adjustable straps)</span>	101
JIV040004	<input type="checkbox"/>	Deep contour - active <span style="float: right;">(horizontal tension adjustable straps)</span>	133
JIV040005	<input type="checkbox"/>	Deep contour - passive <span style="float: right;">(horizontal tension adjustable straps)</span>	133
JIV040009	<input type="checkbox"/>	Omit back upholstery <span style="float: right;">(Chair supplied with NO back upholstery option)</span>	-54



Shallow Contour  
JIV040001 /  
JIV040002 /  
JIV040003



Deep Contour Active  
JIV040004












Deep Contour Passive  
JIV040005

**COMFORT BACK CUSHION WIDTH**








JIV040101	<input type="radio"/>	40cm (15.5")	JIV040106	<input type="radio"/>	42cm (16.5")	JIV040102	<input type="radio"/>	44cm (17")
JIV040107	<input type="radio"/>	46cm (18")	JIV040103	<input type="radio"/>	48cm (19")	JIV040108	<input type="radio"/>	50cm (20")
JIV040104	<input type="radio"/>	52cm (20.5")	JIV040109	<input type="radio"/>	54cm (21")	JIV040105	<input type="radio"/>	56cm (22")

**LATERAL SUPPORTS**

JIV030601	<input type="checkbox"/>	Lateral supports short (pair)	242
JIV030602	<input type="checkbox"/>	Lateral supports long (pair)	242

							£
<b>THIRD PARTY BACKRESTS &amp; INTERFACES</b>							
JIV010302	<input type="checkbox"/>	Jay J3 two point fixation kit					101
See Pages 8 and 9 to select J3 Backrest and options. Additional options can be found on the J3 prescription form							
<b>JAY J3 BACKRESTS AND OPTIONS</b>							
See Pages 8 and 9 to select J3 Backrest and options. Additional options can be found on the J3 prescription form							
<b>POSITIONING STRAPS</b>							
JIV020901	<input type="radio"/>	5cm (2") aircraft buckle					
JIV090056	<input type="radio"/>	5cm (2") aircraft buckle - 190cm (Long)					
JIV020902	<input type="checkbox"/>	5cm (2") padded aircraft buckle					47
<b>HEADRESTS</b>							
JIV040501	<input type="radio"/>	Curved headrest	20cm x 11cm (8" x 4.5")				
JIV040502	<input type="radio"/>	Curved headrest	23cm x 11cm (9" x 4.5")				
JIV040503	<input type="checkbox"/>	Contoured (winged) headrest					79
JIV040504	<input type="checkbox"/>	Jay Comfort: standard comfort seating headrest					91
JIV040505	<input type="checkbox"/>	Jay Comfort: flexible comfort seating headrest (winged style)					111
	Curved Headrest JIV040501 / JIV040502		Contoured Headrest JIV040503		JAY Comfort: Std Comfort headrest JIV040504		JAY Comfort: Winged headrest JIV040505
<b>BACKREST RECLINE OPTIONS</b>							
JIV040601	<input type="radio"/>	Manual adjustable backrest	(-5° to 12°, preset angle at 5°)				
JIV040602	<input type="checkbox"/>	Manual recline gas spring system	(-10° to 30°, require push handle option)				207
JIV040603	<input type="checkbox"/>	Power recline	(-10° to 40°) - Not available with VR2 control system				372
JIV040604	<input type="checkbox"/>	Power recline c/w mechanical shear reduction	(-10° to 40°) - Not available with VR2 control system				690
JIV090503	<input type="checkbox"/>	Height adjustable push handles					91
<b>ARMREST TYPE</b>							
JIV050001	<input type="radio"/>	Single post armrest	(skirt guard to be selected)				
JIV050002	<input type="checkbox"/>	Reclining and flip-back armrest - High	(width and 32 - 38 cm height adjustable)				313
JIV040141	<input type="checkbox"/>	Reclining and flip-back armrest - Low	(width and 27 - 32 cm height adjustable)				313
<b>SKIRT GUARDS</b> <span style="float: right;">(not available with reclining armrest)</span>							
JIV050102	<input type="radio"/>	Long	36cm (14")	(plexi glass)			
JIV050103	<input type="checkbox"/>	Jay Comfort padded skirt guard	36cm (14")				47
<b>ARMPADS - LEFT AND RIGHT HAND SIDE</b>							
 Standard Full  Standard Comfort  JAY Comfort Padded  JAY Comfort P-Shaped  Hemi Armpads	<b>Standard armpads</b>		<b>Left Right</b>				
	Full length	36cm (14")	JIV050201	<input type="radio"/>	<input type="radio"/>	JIV050202	
	Comfort armpads	41cm (16")	JIV050203	<input type="radio"/>	<input type="radio"/>	JIV050204	
	<b>JAY COMFORT armpads</b>						
	Padded (vinyl) upholstery	full length 35cm (13.5")	JIV050251	<input type="checkbox"/>	<input type="checkbox"/>	JIV050252	47 each
	P-Shaped skai upholstery	full length 35cm (13.5")	JIV050253	<input type="checkbox"/>	<input type="checkbox"/>	JIV050254	47 each
	<b>Hemi Armpads</b> <span style="float: right;">(Hemi Armpad &amp; hand pad must be chosen together)</span>						
	Hemi forearm pad	small 23cm (9")	JIV050205	<input type="checkbox"/>	<input type="checkbox"/>	JIV050206	79 each
	Hemi forearm pad	medium 28cm (11")	JIV050207	<input type="checkbox"/>	<input type="checkbox"/>	JIV050208	79 each
	Hemi forearm pad	large 33cm (13")	JIV050209	<input type="checkbox"/>	<input type="checkbox"/>	JIV050210	79 each
Flat hand pads	small 11 x 18cm (4.5" x 7")	JIV050281	<input type="checkbox"/>	<input type="checkbox"/>	JIV050282	79 each	
Flat hand pads	large 14 x 23cm (5.5" x 9")	JIV050283	<input type="checkbox"/>	<input type="checkbox"/>	JIV050284	79 each	
<b>CENTRE MOUNT LEGRESTS</b> <span style="float: right;">(Inc footboard)</span>							
JIV060001	<input type="radio"/>	Fixed centre mount legrest	Short	20.4 cm - 31.6 cm ( 8" - 12.5")			
JIV060002	<input type="radio"/>	Fixed centre mount legrest	Medium	37 cm - 47 cm (13.5" - 18.5")			
JIV060003	<input type="radio"/>	Fixed centre mount legrest	Long	42 cm - 55 cm (16.5" - 22")			
JIV060004	<input type="checkbox"/>	Powered centre mount legrest	(37 - 49cm fixed mechanical adjustment + 14cm via motor:)				982

			£
<b>SWING AWAY LEGRESTS</b>			<i>(width adjustable legrest mount; 8 cm (3") each side)</i>
JIV060101	<input type="radio"/> 70° s/a legrest		
JIV060102	<input type="checkbox"/> 75° s/a legrest tapered	<i>( Only with footboard, reduces inner legrest width by 2")</i>	101
JIV060104	<input type="checkbox"/> Manual elevating adult	<i>(c/w standard foam calf pads)</i>	239
JIV060105	<input type="checkbox"/> Powered elevating	<i>(c/w standard foam calf pads)</i>	451
JIV060106	<input type="checkbox"/> Powered articulating	<i>(c/w standard foam calf pads)</i>	621
JIV060107	<input type="checkbox"/> Auto connectors for powered s/a legrests		101
JIV060109	<input type="checkbox"/> Skai-vinyl calf pads		101
<b>Extension tubes for swing-away legrests:</b>			<i>(per pair)</i>
JIV060201	<input type="radio"/> Short 33 - 41cm (13" - 16")		
JIV060202	<input type="radio"/> Medium 39 - 46cm (15" - 18")		
JIV060203	<input type="radio"/> Long 43 - 53cm (17" - 21")		
<b>FOOTPLATES FOR SWING-AWAY LEGRESTS</b>			<i>(per pair)</i>
JIV060302	<input type="radio"/> Aluminium angle adjustable		
JIV060303	<input type="checkbox"/> Foot board <i>(angle adjustable)</i>	<i>(not available on 90° legrest)</i>	106
JIV060401	<input type="checkbox"/> Extra legrest strap		22
<b>VR2 CONTROLS</b>			<i>(Not Available with 13Kph Motors - Rnet must be selected to deliver increased Power for Kerb/Hill Climbing)</i>
JIV100002	<input type="radio"/> VR2 Joystick: with actuator and light control	<i>(lights &amp; indicators plus max. 2 actuators)</i>	
<b>Dual/Attendant control options</b>			
JIV100201	<input type="checkbox"/> VR2 dual control		398
JIV100202	<input type="checkbox"/> VR2 JSM attendant only		101
	VR2 User Joystick JIV010002		VR2 Dual Control Joystick JIV010201
<b>R-NET 120A CONTROLS</b>			<i>(Required for over 2 actuators &amp; users requiring max torque &amp; performance)</i>
JIV100102	<input type="checkbox"/> R-net control system	<i>(Must be selected if R-Net control system required)</i>	525
JIV100152	<input type="radio"/> R-net LED joystick	<i>(Incl. L&amp;I and actuator control option)</i>	
JIV100153	<input type="checkbox"/> R-net LCD colour joystick	<i>(Integrated colour display with 2 jack inputs for buddy buttons)</i>	59
JIV100154	<input type="radio"/> R-Net Control System - Omit joystick	<i>(When joystick is not required for clients using special controls)</i>	
<b>Dual/Attendant control options</b>			
JIV100201	<input type="checkbox"/> Dual control	<i>(Separate dual control module mounted at rear - Always RH)</i>	398
JIV100202	<input type="checkbox"/> R-net attendant only joystick		101
	R-net LED User Joystick JIV100152		R-net LCD User Joystick JIV100153
	R-net Dual Control Joystick JIV100201		R-net Attendant Only Joystick JIV100202
<b>JOYSTICK MOUNTINGS</b>			
<b>Front control mountings</b>			<i>(for single Joystick selection mount)</i>
JIV100251	<input type="radio"/> Swing-away		
JIV100252	<input type="checkbox"/> Centre bar control mount		153
JIV100253	<input type="checkbox"/> Swing-away tray mounted centre control - R-Net LED	<i>(not available on 22" width)</i>	525
JIV110143	<input type="checkbox"/> Swing-away tray mounted centre control - R-Net LCD	<i>(not available on 22" width)</i>	1460
Joystick mounting side			
JIV100271	Left <input type="radio"/>	Right <input type="radio"/>	JIV100272
<b>LIGHTS &amp; INDICATORS</b>			
JIV110201	<input type="radio"/> LED (StVZO approved)		

								£	
<b>R-NET OMNI DISPLAY &amp; SPECIAL CONTROL INPUT MODULE</b>								(with IR functionality)	
JIV100401	<input type="checkbox"/>	R-net OMNI colour display input module					(includes Special Control Interface)	796	
JIV100550	<input type="checkbox"/>	OMNI in scanning mode						796	
JIV100422	<input checked="" type="checkbox"/>	R-NET OMNI Flexi mount only					(30 cm (12") long flexible height extension)		
		JIV100431	Left	<input type="checkbox"/>	<input type="checkbox"/>	Right	JIV100432		
<b>R-NET OMNI ACCESSORIES</b>									
JIV100503	<input type="checkbox"/>	R-net Bluetooth module					(Compatible with Windows PC and Android systems only)	313	
B4M	<input type="checkbox"/>	R-net Input/Output module					(Build position selected as B4M process)	525	
JIV090516	<input type="checkbox"/>	Bluetooth wireless mouse receiver						59	
	<b>R-net OMNI + Flexi Mount JIV100401</b>  Provides Visual Display, infrared commands (Controls TV, DVD, Stereo, etc) & Special Control interface						<b>R-net Bluetooth Module JIV100503</b>  Allows users to control Windows based PC's, Tablets and Phones from joystick (mounted in base)		
<b>R-NET SPECIAL CONTROLS</b>									
(R-net OMNI required from previous page; 2 buddy buttons recom.)									
JIV100551	<input type="checkbox"/>	Switch-It Versaguide Joystick as hand control						1539	
JIV100552	<input type="checkbox"/>	Switch-It Micro Guide as hand control						2069	
JIV100553	<input type="checkbox"/>	Switch-It Micro Pilot (0 way) joystick as hand control						2281	
JIV100602	<input type="checkbox"/>	Swing-away front mount							
JIV100651	Left	<input type="checkbox"/>	Right	<input type="checkbox"/>			JIV100652		
	<b>R-Net OMNI in Scanning Mode JIV100550</b>			<b>Proportional Versaguide Joystick JIV100551</b>			<b>Switch-It Micro Guide Joystick JIV100552</b>		<b>Switch-It Micro Pilot 0-Way Joystick JIV100553</b>
<b>R-NET HEAD RELATED INPUT DEVICES</b>									
(if you have already chosen a standard headrest then Switch- It head array is not possible)									
JIV100701	<input type="checkbox"/>	Switch-It 3 switch head array / adult			(5 switch and other versions available via B4M)	(OMNI required, 2 BB recommended)	3130		
JIV100703	<input type="checkbox"/>	Switch-It Micro Pilot (0 way) joystick as chin control					(OMNI required, 2 BB recommended)	2281	
JIV100704	<input type="checkbox"/>	Proportional chin control			(proportional joystick; OMNI required; 2 BB recommended)		796		
JIV100705	<input type="checkbox"/>	Sip & Puff control			(tube and mouthpiece; OMNI required; 2 BB recommended)		266		
JIV100706	<input type="checkbox"/>	Switch-It Micro Guide as chin control			(OMNI required, 2 BB recommended)		2069		
JIV100707	<input type="checkbox"/>	Switch-It proportional head control			(OMNI required, 2 BB recommended)		2918		
JIV100708	<input type="checkbox"/>	Switch-It Versaguide Joystick as chin control			(OMNI required, Switch-It flex switch kit incl.)		1645		
	<b>Switch-It Head Array JIV100701</b>			<b>Proportional Versaguide Joystick as chin control JIV100708</b>					
<b>CHIN CONTROL SWING AWAY OPTIONS</b>									
JIV100801	<input type="checkbox"/>	Manual s/a option						207	
	<input type="checkbox"/>	Powered s/a option					(choose your setting below)		
JIV100802	<input type="checkbox"/>	Operated with buddy button			(operated via one buddy button - select colour below)		1539		
JIV100803	<input type="checkbox"/>	Operated with ribbon switch						1645	
	Left	Right	Swing-away direction						
JIV100851	<input type="checkbox"/>	<input type="checkbox"/>	JIV100852						
<b>Buddy button switches</b>									
JIV100901	<input type="checkbox"/>	Red	x 1				79 each		
<b>Piko button switches</b>									
JIV100909	<input type="checkbox"/>	Red	x 1	JIV100913	<input type="checkbox"/>	Red	x 2	79 each	
JIV100910	<input type="checkbox"/>	Green	x 1	JIV100914	<input type="checkbox"/>	Green	x 2	79 each	
JIV100911	<input type="checkbox"/>	Yellow	x 1	JIV100915	<input type="checkbox"/>	Yellow	x 2	79 each	
JIV100912	<input type="checkbox"/>	Blue	x 1	JIV100916	<input type="checkbox"/>	Blue	x 2	79 each	

● = Standard    ○ = No cost option    □ = Cost option			£
<b>REAR or FRONT WHEEL DRIVE</b>			
<b>CASTOR WHEELS</b> <i>(delivered always as a pair)</i> <span style="float:right"><i>(JIVE R/F only)</i></span>			
JIV080001	○	9" x 2.5" Solid	
JIV080002	○	9" x 2.5" Pneumatic	
JIV080003	□	10" x 3" Solid <i>(Not available with 13 kph motors and 6 kph FWD / RWD. Not recommended on RWD base)</i>	79
JIV080004	□	10" x 3" Pneumatic <i>(Not recommended on RWD base)</i>	79
<b>MID WHEEL DRIVE</b>			
JIV080101	●	6" solid castor wheel on front and rear	
<b>KERB CLIMBER</b>			
JIV080501	□	Kerb climber <i>(Jive R2 only)</i>	153
<b>DRIVE BASE</b> <i>(Max kerb climbing is affected by many factors incl user weight, speed selection etc. Contact Tech Serv Ctr for more info)</i>			
<b>MOTORS/SPEED</b>			
JIV090002	○	6kph (4mph) (4-pole) high torque motors for optimal kerb climbing	
JIV090003	○	10kph (6mph) (4-pole) high torque motors for optimal kerb climbing	
JIV090004	□	13kph (8mph) (4-pole) <i>(Not available on JIVE F)</i>	313 <i>(Not available with VR2 - Rnet must be selected)</i>
<b>DRIVE WHEELS</b>			
JIV090101	○	14" Pneumatic "V-type" <i>(Recommended for better kerb climbing)</i>	
JIV090102	○	14" Solid "V-type"	
JIV090103	○	14" Pneumatic "knobby" <i>(Recommended for better kerb climbing)</i>	
JIV090104	○	14" Solid "Knobby"	
<b>BATTERIES</b>			
JIV090201	○	60 Ah maintenance free	
JIV090202	□	70 Ah maintenance free <i>(Recommended for 13Kph Motors where battery power is drawn faster)</i>	101
<b>CHARGER</b>			
JIV090301	○	220V/ 8A automatic charger (UK)	
JIV090303	□	220V/ 10A automatic charger (UK) <i>(Must be ordered with 13kph speed - charges batteries faster)</i>	153
<b>JIVE F/ JIVE R2</b>		<b>JIVE MWD</b>	
<b>(body panel colour, frame is always black)</b>		<b>(body panel colour, frame is always silver)</b>	
JIV010101	○	White RAL 9010	JIV010101 ○ White RAL 9010
JIV010102	○	Red	JIV010102 ○ Red
JIV010103	○	Blue	JIV010103 ○ Blue

		£
<b>ACCESSORIES</b>		
JIV090501	<input type="checkbox"/> Crutch holder	27
JIV090502	<input type="radio"/> Rear view mirror <i>(Always required with 10Kph Chairs)</i>	
JIV090512	<input type="checkbox"/> Swing-away tray table <i>(available from 36 cm (14") to 51 cm (20") seat width)</i>	153
JIV090504	<input type="checkbox"/> 24/12V DC/ DC converter	207

<b>JOYSTICK KNOBS</b>		
JIV090505	<input type="checkbox"/> Golf ball joystick knob	27
JIV090506	<input type="checkbox"/> T bar B38 joystick knob	27
JIV090507	<input type="checkbox"/> Soft ball joystick knob	27
JIV090508	<input type="checkbox"/> Ball joystick knob	27
JIV090509	<input type="checkbox"/> Chin cup joystick knob	27
JIV090510	<input type="checkbox"/> Stick joystick knob	27
JIV090511	<input type="checkbox"/> Mushroom joystick knob	27

<b>TRANSPORTATION OPTIONS</b>		
JIV090539	<input checked="" type="checkbox"/> Prepared for Dahl Docking Station option <i>(Jive F and M via B4M only. Jive R2 docking preparation as standard)</i>	

**For special requirements, please complete Built-4-Me prescription form**



<b>JIVE - base price</b>	£
<b>Option Price</b>	£
<b>Total Price</b>	£



# J3 Back

### J3 Backrest Range

Complete backrest package for clients up to 136kg  
 Universal one point fixation for rapid application  
 Easily adjustable for assessment & changing needs  
 Wide range of positioning and transfer options  
 Backs and options fully transit tested  
 1 Large & 1 small lumbar included in MT, UT & SH box  
 1 small lumbar in LT box      2 year guarantee



- = 5 days from the EDC
- ◇ = 15 days from the EDC
- = upcharge option


Max. user weight: **136kg**  
 USER WEIGHT  kg  
 Mandatory information to ensure correct configuration

#### Applications:

For all clients requiring postural back support from mild lateral stability to extreme lateral and posterior stability.  
 Full range of postural back shapes to cover clinical back support at the appropriate backrest height for the user.  
 Adjustable in depth, height, angle and chair width, ensuring rapid and comprehensive client positioning and spine alignment.

### CLIENT MEASUREMENT

1. Assess Lateral Stability Requirement (SC, MC, DC or PDC)
2. Assess Posterior Stability Requirement (MT, UT or SH)
3. Assess appropriate backrest height (Short, Medium, Tall)
4. Determine client trunk width across the sternum; and cross check with backpost width.
5. Assess for options and 'Velcro' positioning inserts



**Item Nr. construction:**  
**J3 SC LT S 41** (This is only an example)  
 J3: J3 Back  
 SC: Shallow Contour  
 LT: Lower Thoracic  
 S: Short (Height cm)  
 41: Chair Backrest Width (cm)

#### Width


Chair backrest widths (Measured from outside of backrest to outside of backrest)

Jay 3 Backrest Width																	
41			46			51			41			46			51		
41 - 46	46 - 51	51 - 56	41 - 46	46 - 51	51 - 56	41 - 46	46 - 51	51 - 56	41 - 46	46 - 51	51 - 56	41 - 46	46 - 51	51 - 56			

#### Back height



Client max trunk width (cm)		£520			£520			£570			£570		
Upper Thoracic (UT)	Short (S) - 42 Med (M) - 46 Tall (T) - 50	38	43	48	34,5	39,5 (1)	44,5	33	38 (1)	43	33	38	43
		○	○	○	○	○	○	○	○	○	○	○	○



Client max trunk width (cm)		£520			£520			£570			£570		
Shoulder Height (SH)	Short (S) - 53 Med (M) - 57 Tall (T) - 61	38	43	48	34,5	39,5 (1)	44,5	33	38 (1)	43	33	38	43
		○	○	○	○	○	○	○	○	○	○	○	○

(1) Not available with laterals

#### Conversion Table

cm	17.0	20.0	24.0	30.0	34.0	36.0	38.0	41.0	42.0	44.0	46.0	51.0	53.0	56.0	61.0
inches	6.5	8.0	9.5	12.0	13.5	14.0	15.0	16.0	16.5	17.5	18.0	20.0	21.0	22.0	24.0



JAY© 3 BACK ACCESSORIES

= upcharge option

ACCESSORIES				Unit Price
JAY090001	<input type="checkbox"/>	Backrest depth extension		£35
JAY090002	<input type="checkbox"/>	4 Point hardware kit	<b>NEW!</b>	£147

VANITY FLAP	41 cm	46 cm	51 cm	Unit Price
Vanity flap (covers gap between cushion and backrest)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£47
	JAY090003	JAY090004	JAY090005	

**HEADREST** (only available where headrest is indicated on the first page) **Crash tested** Unit Price

Please choose 1 mount and 1 pad to have a complete headrest with pad

**Mount**

JAY090006	<input type="checkbox"/>	Jay standard headrest (bracket & mount only)	£80
JAY090007	<input type="checkbox"/>	Jay flip back headrest (bracket & mount only)	£180

**Pad**

JAY090008	<input type="checkbox"/>	Jay standard pad	£37
JAY090009	<input type="checkbox"/>	Jay contour pad	£37
JAY090010	<input type="checkbox"/>	Universal headrest mount (required to use any other headrest not listed in this order form)	£37

**LATERAL SUPPORT** (only available where laterals are indicated on the first page) Unit Price

Left Hand lateral		Depth					
cm		10	12	15			
Height	10	<input type="checkbox"/> JAY090011	<input type="checkbox"/> JAY090012	<input type="checkbox"/> JAY090013	←	Fix lateral JAY090017	£42
	15	<input type="checkbox"/> JAY090014	<input type="checkbox"/> JAY090015	<input type="checkbox"/> JAY090016	←	Fix lateral JAY090018	£42
		<input type="checkbox"/> JAY090019	←	Swing-away option		£79	
Right Hand lateral		Depth					
cm		10	12	15			
Height	10	<input type="checkbox"/> JAY090020	<input type="checkbox"/> JAY090021	<input type="checkbox"/> JAY090022	←	Fix lateral JAY090026	£42
	15	<input type="checkbox"/> JAY090023	<input type="checkbox"/> JAY090024	<input type="checkbox"/> JAY090025	←	Fix lateral JAY090027	£42
		<input type="checkbox"/> JAY090028	←	Swing-away option		£79	

When ordering 2 laterals Left and Right, you need to pick the sizes in both matrices above.

**Total Price for laterals**

If a swing-away lateral is required, first select the height and the depth, then tick the swing-away option box below the matrix

POSITIONING KIT		CM width	41 cm	46 cm	51 cm	Unit Price
Spine Align Mini Fit Kit I			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£28
			JAY090029	JAY090030	JAY090031	
Spine Align Mini Fit Kit II			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£43
			JAY090032	JAY090033	JAY090034	
JAY090035	<input type="checkbox"/>	Jay 3 Back Positioning Box				£137
JAY090036	<input type="checkbox"/>	Jay 3 Back Hardware Box				£431

CHEST STRAP		(only available where harnesses are indicated on the first page)			Medium	Large	Extra large	Unit Price
Padded anterior strap			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			£143
			JAY090037	JAY090038	JAY090039			
Centre open harness			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			£143
			JAY090040	JAY090041	JAY090042			

SPARES		(separate to product - separate box or delivery)			CM width	41 cm	46 cm	51 cm	Unit Price
Small lumbar		<i>(1 unit included as std with product LT, MT, UT &amp; SH)</i>				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£8.00
			JAY090043	JAY090044	JAY090045				
Large lumbar support		<i>(1 unit included as std with product MT, UT &amp; SH)</i>				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£8.00
			JAY090046	JAY090047	JAY090048				
Wedge			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				£8.00
			JAY090049	JAY090050	JAY090051				

		Unit price	Quantity	Total Price
JAY090052	<input type="checkbox"/>	Curved shim	£8.00	
JAY090053	<input type="checkbox"/>	Flat shim	£8.00	
JAY090054	<input type="checkbox"/>	Posterior/lateral pelvic	£11.00	
JAY090055	<input type="checkbox"/>	Lateral shim	£4.00	
JAY090056	<input type="checkbox"/>	Small bean bag	£8.00	
JAY090057	<input type="checkbox"/>	Large bean bag	£8.00	