

Invoice to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Order date:	
Purchase order:	

Deliver to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Marked for:	

Base Pricing
£1,995.00

E-mail: orders@sunmed.co.uk

Spare Parts: <http://parts.sunrisemedical.eu>

Fax: +44 (0)845 605 66 89

Maximum user weight limit 136 Kg (21.5 stone)

Part No 21400000-000



Standard features	<ul style="list-style-type: none"> ● Width, height & depth adjustable seat ● Flip back, angle, height & depth adjustable arms ● Height and angle adjustable backrest ● Tension adjustable backrest sling 	<ul style="list-style-type: none"> ● Fixed centre mount legrest ● Adjustable lap strap ● Transportation approved
<ul style="list-style-type: none"> ● Mid wheel drive base ● 10" solid drive wheels ● 36Ah batteries 		

STANDARD SEAT: FRAME / SEAT SIZE

Width Setting		Depth setting	
HULSW41	<input type="radio"/> 41cm (16")	HULSD41	<input type="radio"/> 41cm (16")
HULSW43	<input type="radio"/> 43cm (17")	HULSD43	<input type="radio"/> 43cm (17")
HULSW46	<input type="radio"/> 46cm (18")	HULSD46	<input type="radio"/> 46cm (18")
HULSW48	<input type="radio"/> 48cm (19")	HULSD48	<input type="radio"/> 48cm (19")
HULSW51	<input type="radio"/> 51cm (20")	HULSD51	<input type="radio"/> 51cm (20")

Seat frame is fully adjustable from 41cm to 51cm (16" - 20") in width and depth. Independently adjustable in 2.5cm (1") increments

● = Standard ○ = No cost option □ = Cost option

SEAT & CUSHION

HUL020000	● Firm seat board	
HUL020201	<input type="radio"/> 5cm (2") standard foam cushion	(41cm or 48cm option sized to seat)
JCL020003	<input type="checkbox"/> JAY COMFORT - Standard cushion - stretched rubber upholstery	(thickness 8cm) 137
JCL020007	<input type="checkbox"/> JAY COMFORT - Optimum comfort - stretched rubber upholstery	(anatomically formed - thickness 9cm) 176
HUL020203	<input type="radio"/> Omit cushion	

SEAT TILT / ANGLE

HUL020204	● Manual tilt (3° preset tilt)	(0°, 3° and 6° adjustable)
-----------	--------------------------------	----------------------------

SEAT HEIGHT

HULSH45	<input type="radio"/> 45cm (18") seat to floor height	(with 3° seat angle/tilt)	(0° seat tilt: 41.5cm / 6° seat tilt 48.5cm)
HULSH47	<input type="radio"/> 47.5cm (19") seat to floor height	(with 3° seat angle/tilt)	(0° seat tilt: 44cm / 6° seat tilt 51cm)

BACKREST

HUL030000	● Mechanical fixed recline back	(-4° to 12°)
HUL030001	<input type="radio"/> Tension adj. sling backrest (48cm /19" height)	(2.5cm up and down adjustment, with standard foam cushion only)
HUL030009	<input type="radio"/> Omit backrest upholstery	

JAY COMFORT BACKRESTS

JCL030004	<input type="checkbox"/> Shallow contour (standard velour) - 50cm back height	(hard back)	135
JCL030008	<input type="checkbox"/> Shallow contour (standard velour) - 57cm back height	(vertical elastic straps)	196
JCL030006	<input type="checkbox"/> Shallow contour (standard velour) - 57cm back height	(horizontal tension adjustable straps)	196
JCL030012	<input type="checkbox"/> Deep contour - passive - 57cm back height	(horizontal tension adjustable straps) (4" integrated lateral support)	227
JCL030010	<input type="checkbox"/> Deep contour - active - 57cm back height	(horizontal tension adjustable straps) (4" integrated lumbar lateral support)	227

● = Standard ○ = No cost option □ = Cost option

				£	
LATERAL SUPPORT: PAIR				(Width, depth, angle height adjustable)	
HUL030007	<input type="checkbox"/>	Lateral supports: Short	(Only available with tension adj backrests)	(5" Depth from backpost to centre of pad)	283
HUL030008	<input type="checkbox"/>	Lateral supports: Long	(recommended with JAY Comfort backrests)	(10" Depth from backpost to centre of pad)	283
HEADREST					
HUL090001	<input type="checkbox"/>	Comfort headrest (winged style with adj. side pads)		(Mounted via push handles)	153
HUL090002	<input type="checkbox"/>	JAY Comfort: standard comfort seating headrest		(Mounted via Comfort backrest - only available with Comfort back)	173
ARMREST					
	●	Flip-up armrest (height, angle adjustable)			
LEGRESTS & FOOTPLATES					
Fixed centre mount legrests & flip-up footboard					
HUL050000	○	Fixed centre mount legrest		(+/- 1.5cm height adjustable)	
Swing-away legrests					
HUL050001	<input type="checkbox"/>	80 degree swing-away legrests		(34 - 50cm Lower leg length)	98
HUL050002	<input type="checkbox"/>	Manual elevating legrests			180
Fixed angle footplates					
HUL050006	○	Fixed angle footplates			
HUL050003	<input type="checkbox"/>	Angle adjustable footplates			52
Amputee Support					
HUL050004	<input type="checkbox"/>	Amputee support - RH (Uni style)			95
HUL050005	<input type="checkbox"/>	Amputee support - LH (Uni style)			95
VR2 CONTROLS					
HUL110000	●	VR2 remote			
	●	Remote control mounting	HUL110002	RH ○ ○ LH HUL110003	
MOTOR SPEED					
HUL100000	○	3.2kph (2mph indoor setting)		(6kph motor with reduced speed via programming)	
HUL100001	○	6kph (4mph)			
DRIVE BASE FEATURES & ACCESSORIES					
	●	10" Solid drive wheels	6" solid front and rear castors		
	●	Red Shroud			
	●	Transportation tie down points			
	●	Lap belt			
	●	36Ah batteries			
HUL030200	●	5A charger			
HUL090005	<input type="checkbox"/>	Swing-away tray table			153
Quickie HULA MWD Base Price				£1,995.00	
Option Price				£	
VAT				£	
Total Price				£	