

Invoice to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Order date:	
Purchase order:	

Deliver to:	
Name:	
Street:	
Town:	
Country:	
Tel:	
Marked for:	

Standard Features		
• Power double tilt module	• Chest fixation (2-point strap)	• Swing-away control mount
• 0° - 22° Posterior Power tilt	• Single post armrest	• LED lights and indicators
• 0° - 80° Anterior Power tilt	• Standard full length armpads	• R-Net complex seating / lighting
• Powered lift 20cm	• R-Net 120A Controller	• 6kph (4mph) 4 pole motors
• 0° to 85° powered recline with 100/110/120 mm anti-shear mechanism	• R-Net LCD Joystick	• Pneumatic tyres with V tread
• Powered centre mount legrest (Incl. Knee fixation)		• 60Ah batteries
		• 8A Charger

Fax: +44 (0)845 605 66 89

E-mail: orders@sunmed.co.uk

Spare Parts: <http://parts.sunrsemmedical.eu>

JIVE UP: £14,950

Max. user weight: 120kg (19 stone)

SEAT DIMENSION JAY COMFORT

Seat frame width is 41 cm. The Seat width can be adjusted between 40 and 52 cm via the armrests, based on the seat/back cushion configuration. The Seat depth can be adjusted between 41 and 55 cm via backrest movement

Seat dimensions also reflect aluminium seat plate option. Seat plate option to be selected at the Jay Comfort Cushion section (JIV020109).

Cushion width

cm	40	42	44	46	48	50		
Seat width	40	○ JUP020001	○ JUP020008	○ JUP020015	○ JUP020022	○ JUP020029	○ JUP020036	Seat cushion depth 48 cm
	42	○ JUP020002	○ JUP020009	○ JUP020016	○ JUP020023	○ JUP020030	○ JUP020037	
	44	○ JUP020003	○ JUP020010	○ JUP020017	○ JUP020024	○ JUP020031	○ JUP020038	
	46	○ JUP020004	○ JUP020011	○ JUP020018	○ JUP020025	○ JUP020032	○ JUP020039	
	48	○ JUP020005	○ JUP020012	○ JUP020019	○ JUP020026	○ JUP020033	○ JUP020040	Seat cushion depth 55 cm
	50	○ JUP020006	○ JUP020013	○ JUP020020	○ JUP020027	○ JUP020034	○ JUP020041	
	52	○ JUP020007	○ JUP020014	○ JUP020021	○ JUP020028	○ JUP020035	○ JUP020042	

SEAT HEIGHT

Seat height

cm	45	48
Jive M	○ SH45	○ SH48

Seat height measured floor to seat pan at 0° tilt and without seat cushion

SEAT TO BASE INTERFACE

JUP130010	Power stand-up modules	
	● Power double tilt module	(counts as 2 actuators)
	● 0° - 22° Posterior Power tilt	
	● 0° - 80° Anterior Power tilt	
	● Powered lift 20cm	
	● 0° to 85° powered recline with 100/110/120 mm anti-shear mechanism	
	● Powered centre mount legrest (incl. Knee fixation)	

JAY CUSHION SIZE

SEAT WIDTH

		Width Coding	JL 40	JL 42	JL 44	JL 46	JL 48	JL 50
Depth Coding	cm	40	42	44	46	48	50	
JL __ 40	40	○	○	○	○(1)	○(1)		
JL __ 42	42	○	○	○	○	○		
JL __ 44	44	○	○	○	○	○		
JL __ 46	46	○	○	○	○	○	○	
JL __ 50	50				○	○	○	

Size Configuration - COMPLETE PRODUCT

Conversion Table

cm	40.0	42.0	44.0	46.0	48.0	50.0	56.0
inches	15.6	16.4	17.2	18.0	18.8	19.5	21.9



(1) Not available with Soft Combi P or Easy Fluid Cushion

JAY CUSHION OPTIONS

JAY020004	<input type="checkbox"/>	Jay Lite Standard version	(Micro climatic cover, Aluminium seat plate to be selected)	195
JAY020005	<input type="checkbox"/>	Jay Easy Fluid with Flat Base	(Incontinence cover, Aluminium seat plate to be selected)	280
JAY020001	<input type="checkbox"/>	Jay Easy Visco with Flat Base	(Incontinence cover, Aluminium seat plate to be selected)	120
JAY020003	<input type="checkbox"/>	Jay Soft Combi P	(Incontinence cover, Aluminium seat plate to be selected)	70
JUP020109	<input type="radio"/>	Seat Plate Aluminium	(supplied in the size of the cushion chosen)	

JAY COMFORT SEAT CUSHION OPTIONS

JCL020003	<input type="checkbox"/>	Standard cushion - stretched rubber upholstery	(thickness 8cm)	141
JCL020007	<input type="checkbox"/>	Optimum comfort - stretched rubber upholstery	(anatomically formed) (thickness 9cm)	181
JUP020110	<input type="radio"/>	Omit cushion		



JAY COMFORT BACK CUSHION OPTION

(57cm back height - Adj. up to 70cm)

JCL030004	<input type="radio"/>	Shallow contour (standard velour)	(50 cm back height)	(hard back)	139
JCL030008	<input type="checkbox"/>	Shallow contour (standard velour)	(57 cm back height)	(vertical elastic straps)	202
JCL030006	<input type="checkbox"/>	Shallow contour (standard velour)	(57 cm back height)	(horizontal tension adjustable straps)	202
JCL030010	<input type="checkbox"/>	Deep contour - active	(57 cm back height)	(horizontal tension adjustable straps)	234
JCL030012	<input type="checkbox"/>	Deep contour - passive	(57 cm back height)	(horizontal tension adjustable straps)	234
JUP030209	<input type="radio"/>	Omit back upholstery	(Chair supplied with NO back upholstery option)		
JUP030210	<input type="checkbox"/>	Jay Comfort mounting hardware only			35



Shallow Contour
JCL030004 /
JCL030008 /
JCL030006



Deep Contour Active
JCL030010











Deep Contour
Passive
JCL030012

COMFORT BACK CUSHION WIDTH

BW40	<input type="radio"/>	40cm (15.5")	BW42	<input type="radio"/>	42cm (16.5")	BW44	<input type="radio"/>	44cm (17")
BW46	<input type="radio"/>	46cm (18")	BW48	<input type="radio"/>	48cm (19")	BW50	<input type="radio"/>	50cm (20")

POSITIONING STRAPS

JUP090047	<input type="radio"/>	5 cm (2") aircraft buckle				
JUP090048	<input type="checkbox"/>	5 cm (2") padded aircraft buckle				95
JUP090052	●	Chest fixation (2-point strap)	JUP090054	M	L	JUP090055
JUP090053	<input type="checkbox"/>	2nd Chest fixation (2-point strap)	JUP090056	M	L	JUP090057
				○	○	
				○	○	

										£	
HEADRESTS											
JUP030426	<input type="radio"/>	Curved headrest	20cm x 11cm (8" x 4.5")								
JUP030417	<input type="radio"/>	Curved headrest	23cm x 11cm (9" x 4.5")								
JUP030427	<input type="checkbox"/>	Contoured (winged) headrest								95	
JUP030428	<input type="checkbox"/>	Jay Comfort: standard comfort seating headrest								125	
JUP030429	<input type="checkbox"/>	Jay Comfort: flexible comfort seating headrest (winged style)								150	
JUP030415	<input type="radio"/>	Omit headrest									
	Curved Headrest JUP030426 / JUP030417			Contoured Headrest JUP030427			JAY Comfort: Std Comfort headrest JUP030428			JAY Comfort: Winged headrest JUP030429	
BACKREST RECLINE OPTIONS											
JUP030204	<input type="checkbox"/>	Height adjustable push handles								88	
JUP030601	<input type="checkbox"/>	Height adjustable swing-away laterals (short)								235	
ARMREST TYPE											
JUP040122	<input type="radio"/>	Single post armrest								<i>(width and height adjustable, pair)</i>	
JUP040138	<input type="checkbox"/>	Reclining and flip-back armrest - High								<i>(width and height adjustable, pair)</i> 304	
SKIRT GUARDS											
										<i>(not available with reclining armrest)</i>	
JUP040139	<input type="checkbox"/>	Jay Comfort padded skirt guard	36cm (14")								95
ARMPADS - LEFT AND RIGHT HAND SIDE											
   	Standard armpads		Left Right								
	Full length		42 cm (16,5")	JUP040201	<input type="radio"/>	<input type="radio"/>	JUP040202				
	JAY COMFORT armpads										
	Padded (vinyl) upholstery		full length 40 cm (16")	JUP040203	<input type="checkbox"/>	<input type="checkbox"/>	JUP040204	46 each			
	P-Shaped skai upholstery		full length 40 cm (16")	JUP040205	<input type="checkbox"/>	<input type="checkbox"/>	JUP040206	46 each			
	Hemi Armpads										
	<i>(Hemi Armpad & hand pad must be chosen together)</i>										
	Hemi forearm pad		small 23cm (9")	JUP040207	<input type="checkbox"/>	<input type="checkbox"/>	JUP040208	77 each			
	Hemi forearm pad		medium 28cm (11")	JUP040209	<input type="checkbox"/>	<input type="checkbox"/>	JUP040210	77 each			
	Flat hand pads		small 11 x 18cm (4.5" x 7")	JUP040211	<input type="checkbox"/>	<input type="checkbox"/>	JUP040212	77 each			
Flat hand pads		large 14 x 23cm (5.5" x 9")	JUP040213	<input type="checkbox"/>	<input type="checkbox"/>	JUP040214	77 each				
R-NET 120A CONTROLS											
JUP110004	<input checked="" type="radio"/>	R-net 120A control system								<i>(no on board actuator control output)</i>	
JUP110017	<input checked="" type="radio"/>	R-net LCD colour joystick								<i>(Integrated colour display with 2 jack inputs for buddy buttons)</i>	
Dual control options											
JUP110021	<input type="checkbox"/>	Dual control								<i>(Separate dual control module mounted at rear - Always RH)</i> 386	
JOYSTICK MOUNTINGS											
Front control mountings											
<i>(for single Joystick selection mount)</i>											
JUP110051	<input type="radio"/>	Swing-away									
JUP110143	<input type="checkbox"/>	Swing-away tray mounted centre control - R-Net LCD								<i>(not available on 22" width)</i> 1417	
Joystick mounting side											
JUP110053		Left	<input type="radio"/>	<input type="radio"/>	Right	JUP110052					
5 Switch direct actuator - only with colour LCD JSM (JUP110017)											
JUP110171	<input type="checkbox"/>	5-Switch direct actuator box with paddle switches								<i>(4 predefined and 1 programmable seat position)</i> 450	
JUP110172	<input type="checkbox"/>	5-Switch direct actuator box with buttons								<i>(4 predefined and 1 programmable seat position)</i> 450	

LIGHTS & INDICATORS

JIV110201	<input type="radio"/>	LED (StVZO approved)	
-----------	-----------------------	----------------------	--

R-NET OMNI DISPLAY & SPECIAL CONTROL INPUT MODULE

(with IR functionality)

JUP110101	<input type="checkbox"/>	R-net OMNI colour display input module	(includes Special Control Interface)	773
JUP110102	<input type="checkbox"/>	OMNI in scanning mode	(c/w 10 assignable functions through DB9 interface)	773
JUP110110	<input type="radio"/>	R-NET OMNI Flexi mount only	(30 cm (12") long flexible height extension)	
		JUP110121 Left <input type="radio"/>	<input type="radio"/>	Right JUP110122

R-NET OPTIONAL CONTROLS

JUP110033	<input checked="" type="checkbox"/>	R-NET CxSM (seating module with memory function)		
JUP110031	<input checked="" type="checkbox"/>	R-NET ISM-L (lighting module)		
JUP090068	<input type="checkbox"/>	R-net Bluetooth module	(Compatible with Windows PC and Android systems only)	304
JUP090077	<input type="checkbox"/>	R-net 2nd Bluetooth module	(Compatible with Windows PC and Android systems only)	304
JUP090078	<input type="checkbox"/>	R-net Input/Output module	(7 Switched controllable output commands, only available via B4Me)	350
JUP090069	<input type="checkbox"/>	Bluetooth wireless mouse receiver		57

R-NET SPECIAL CONTROLS

(R-net OMNI required from previous page; 2 buddy buttons recom.)

JUP110130	<input type="checkbox"/>	Switch-It Versaguide Joystick as hand control		1494
JUP110131	<input type="checkbox"/>	Switch-It Micro Guide as hand control		2009
JUP110132	<input type="checkbox"/>	Switch-It Micro Pilot (0 way) joystick as hand control		2215
JUP110150	<input type="radio"/>	Fixed front mount		
JUP110151	<input type="radio"/>	Swing-away front mount		
JUP110162	<input type="radio"/>	Left <input type="radio"/>	Right <input type="radio"/>	JUP110161



R-Net OMNI in Scanning Mode
JIV100550



Proportional Versaguide Joystick
JIV100551



Switch-It Micro Guide Joystick
JIV100552



Switch-It Micro Pilot 0-Way Joystick
JIV100553

R-NET HEAD RELATED INPUT DEVICES

(if you have already chosen a standard headrest then Switch- It head array is not possible)

JUP110133	<input type="checkbox"/>	Switch-It 3 switch head array / adult	(5 switch and other versions available via B4M)	(OMNI required, 2 BB recommended)	3039
JUP110138	<input type="checkbox"/>	Switch-It proportional head control		(OMNI required, 2 BB recommended)	2833



Switch-It Head Array
JUP110133

Buddy button switches

JUP110181	<input type="checkbox"/>	Red	x 1		77 each
-----------	--------------------------	-----	-----	--	---------

Piko button switches

JUP110191	<input type="checkbox"/>	Red	x 1	JUP110192	<input type="checkbox"/>	Red	x 2	77 each
JUP110193	<input type="checkbox"/>	Green	x 1	JUP110194	<input type="checkbox"/>	Green	x 2	77 each
JUP110195	<input type="checkbox"/>	Yellow	x 1	JUP110196	<input type="checkbox"/>	Yellow	x 2	77 each
JUP110197	<input type="checkbox"/>	Blue	x 1	JUP110198	<input type="checkbox"/>	Blue	x 2	77 each

SPEED

JIV100710	<input type="radio"/>	6 kph (4mph)	(4 - pole)	
JIV100711	<input type="radio"/>	10 kph (6mph)	(4 - pole)	
JIV100712	<input type="checkbox"/>	13 kph (8mph)	(4 - pole)	(lower torque than 10kph motors) 313

DRIVE WHEELS - JIVE M always with 6" solid castor wheel on front and rear

JIV070032	<input type="radio"/>	14" Pneumatic "V-type"	
JIV070031	<input type="radio"/>	14" Solid "V-type"	
JIV070034	<input type="radio"/>	14" Pneumatic "knobby"	
JIV070033	<input type="radio"/>	14" Solid "Knobby"	

● = Standard ○ = No cost option

□ = Cost option

£

BATTERIES

JIV160001	<input type="radio"/>	60 Ah maintenance free	
JIV160002	<input type="checkbox"/>	70 Ah maintenance free	(Recommended for 13Kph Motors where battery power is drawn faster) 101

CHARGER

JIV160100	<input type="radio"/>	220V/ 8A automatic charger (UK)	
JIV160102	<input type="checkbox"/>	220V/ 10A automatic charger (UK)	(Must be ordered with 13kph speed - charges batteries faster) 153

BODY PANEL COLOUR - frame is always silver /black

JIV140001	<input type="radio"/>	White RAL 9010	
JIV140002	<input type="radio"/>	Red	
JIV140003	<input type="radio"/>	Blue	

ACCESSORIES

JUP090067	<input type="radio"/>	Rear view mirror	(Always required with 10Kph Chairs)
JUP090040	<input type="checkbox"/>	Swing-away tray table	149
JUP090061	<input type="checkbox"/>	24/12V DC/ DC converter	201

JOYSTICK KNOBS

JUP090505	<input type="checkbox"/>	Golf ball joystick knob	26
JUP090506	<input type="checkbox"/>	T bar B38 joystick knob	26
JUP090507	<input type="checkbox"/>	Soft ball joystick knob	26
JUP090508	<input type="checkbox"/>	Ball joystick knob	26
JUP090509	<input type="checkbox"/>	Chin cup joystick knob	26
JUP090510	<input type="checkbox"/>	Stick joystick knob	26
JUP090511	<input type="checkbox"/>	Mushroom joystick knob	26

For special requirements, please complete Built-4-Me prescription form



JIVE - base price	£14,950
Option Price	£
Total Price	£

J3 Back

J3 Backrest Range

Complete backrest package for clients up to 136kg
 Universal one point fixation for rapid application
 Easily adjustable for assessment & changing needs
 Wide range of positioning and transfer options
 Backs and options fully transit tested
 1 Large & 1 small lumbar included in MT, UT & SH box
 1 small lumbar in LT box 2 year guarantee



- = 5 days from the EDC
- ◇ = 15 days from the EDC
- = upcharge option

Max. user weight: **136kg**
 USER WEIGHT kg

Mandatory information to ensure correct configuration

Applications:

For all clients requiring postural back support from mild lateral stability to extreme lateral and posterior stability.
 Full range of postural back shapes to cover clinical back support at the appropriate backrest height for the user.
 Adjustable in depth, height, angle and chair width, ensuring rapid and comprehensive client positioning and spine alignment.

CLIENT MEASUREMENT

1. Assess Lateral Stability Requirement (SC, MC, DC or PDC)
2. Assess Posterior Stability Requirement (MT; UT or SH)
3. Assess appropriate backrest height (Short, Medium, Tall)
4. Determine client trunk width across the sternum; and cross check with backpost width.
5. Assess for options and 'Velcro' positioning inserts

5 cm

5 cm

Shallow Contour / SC
Moderate Lateral Support

8 cm

8 cm

Mid Contour / MC
Medium Lateral Support

15 cm

15 cm

Deep Contour / DC
Full Lateral Support

15 cm

15 cm

Posterior Deep Contour / PDC
PDC Full Lateral support

Item Nr. construction:
J3 SC LT S 41 (This is only an example)
J3: J3 Back
SC: Shallow Contour
LT: Lower Thoracic
S: Short (Height cm)
41: Chair Backrest Width (cm)

Width

Chair backcane widths (Measured from outside of backcane to outside of backcane)

Back height

Jay 3 Backrest Width																	
41			46			51			41			46			51		
	41 - 46	46 - 51	46 - 51	51 - 56	51 - 56	41 - 46	46 - 51	51 - 56	41 - 46	46 - 51	51 - 56	41 - 46	46 - 51	51 - 56	41 - 46	46 - 51	51 - 56

Laterals, headrest & harnesses available													
		£520			£520			£570			£570		
Backrest Height		38	43	48	34,5	39,5 (1)	44,5	33	38 (1)	43	33	38	43
Upper	(S) - 42	○	○	○	○	○	○	○	○	○	○	○	○
Thoracic	(M) - 46	○	○	○	○	○	○	○	○	○	○	○	○
(UT)	(T) - 50	○	○	○	○	○	○	○	○	○	○	○	○

Laterals, headrest & harnesses available													
		£520			£520			£570			£570		
Backrest Height		38	43	48	34,5	39,5 (1)	44,5	33	38 (1)	43	33	38	43
Shoulder	(S) - 53	○	○	○	○	○	○	○	○	○	○	○	○
Height	(M) - 57	○	○	○	○	○	○	○	○	○	○	○	○
(SH)	(T) - 61	○	○	○	○	○	○	○	○	○	○	○	○

(1) Not available with laterals

Conversion Table

cm	17.0	20.0	24.0	30.0	34.0	36.0	38.0	41.0	42.0	44.0	46.0	51.0	53.0	56.0	61.0
inches	6.5	8.0	9.5	12.0	13.5	14.0	15.0	16.0	16.5	17.5	18.0	20.0	21.0	22.0	24.0

JAY© 3 BACK ACCESSORIES

= upcharge option

MOUNTING HARDWARE				Unit Price
JAY090060	<input type="checkbox"/>	Jay 3 Back mounting hardware only		£98

VANITY FLAP	41 cm	46 cm	51 cm	Unit Price
Vanity flap (covers gap between cushion and backrest)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£47
	JAY090003	JAY090004	JAY090005	

HEADREST (only available where headrest is indicated on the first page)	Crash tested	Unit Price
Please choose 1 mount and 1 pad to have a complete headrest with pad		

Mount			
JAY090006	<input type="checkbox"/>	Jay standard headrest (bracket & mount only)	£80
JAY090007	<input type="checkbox"/>	Jay flip back headrest (bracket & mount only)	£180

Pad			
JAY090008	<input type="checkbox"/>	Jay standard pad	£37
JAY090009	<input type="checkbox"/>	Jay contour pad	£37

JAY090010	<input type="checkbox"/>	Universal headrest mount (required to use any other headrest not listed in this order form)	£37
-----------	--------------------------	---	-----

LATERAL SUPPORT (only available where laterals are indicated on the first page)	Unit Price
---	------------

Left Hand lateral

Height	Depth							
	cm	10	12			15		
10	<input type="checkbox"/>	JAY090011	<input type="checkbox"/>	JAY090012	<input type="checkbox"/>	JAY090013	← Fix lateral JAY090017	£42
15	<input type="checkbox"/>	JAY090014	<input type="checkbox"/>	JAY090015	<input type="checkbox"/>	JAY090016	← Fix lateral JAY090018	£42
	<input type="checkbox"/>	JAY090019	← Swing-away option					£79

Right Hand lateral

Height	Depth							
	cm	10	12			15		
10	<input type="checkbox"/>	JAY090020	<input type="checkbox"/>	JAY090021	<input type="checkbox"/>	JAY090022	← Fix lateral JAY090026	£42
15	<input type="checkbox"/>	JAY090023	<input type="checkbox"/>	JAY090024	<input type="checkbox"/>	JAY090025	← Fix lateral JAY090027	£42
	<input type="checkbox"/>	JAY090028	← Swing-away option					£79

When ordering 2 laterals Left and Right, you need to pick the sizes in both matrices above. **Total Price for laterals**

If a swing-away lateral is required, first select the height and the depth, then tick the swing-away option box below the matrix

POSITIONING KIT				CM width	41 cm	46 cm	51 cm	Unit Price
Spine Align Mini Fit Kit I					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£28
					JAY090029	JAY090030	JAY090031	
Spine Align Mini Fit Kit II					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£43
					JAY090032	JAY090033	JAY090034	
JAY090035	<input type="checkbox"/>	Jay 3 Back Positioning Box						£137
JAY090036	<input type="checkbox"/>	Jay 3 Back Hardware Box						£431

CHEST STRAP (only available where harnesses are indicated on the first page)				Medium	Large	Extra large	Unit Price
Padded anterior strap				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£143
				JAY090037	JAY090038	JAY090039	
Centre open harness				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£143
				JAY090040	JAY090041	JAY090042	

SPARES (separate to product - separate box or delivery)				CM width	41 cm	46 cm	51 cm	Unit Price
Small lumbar <small>(1 unit included as std with product LT, MT, UT & SH)</small>					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£8.00
					JAY090043	JAY090044	JAY090045	
Large lumbar support <small>(1 unit included as std with product MT, UT & SH)</small>					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£8.00
					JAY090046	JAY090047	JAY090048	
Wedge					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£8.00
					JAY090049	JAY090050	JAY090051	

		Unit price	Quantity	Total Price
JAY090052	<input type="checkbox"/>	Curved shim		£8.00
JAY090053	<input type="checkbox"/>	Flat shim		£8.00
JAY090054	<input type="checkbox"/>	Posterior/lateral pelvic		£11.00
JAY090055	<input type="checkbox"/>	Lateral shim		£4.00
JAY090056	<input type="checkbox"/>	Small bean bag		£8.00
JAY090057	<input type="checkbox"/>	Large bean bag		£8.00

 <p>JAY090056 & JAY090057</p>	 <p>JAY090006</p>	 <p>JAY090007</p>
 <p>JAY090008</p>	 <p>JAY090009</p>	 <p>JAY090010</p>
 <p>Fixed laterals</p>	 <p>Swing-away laterals</p>	 <p>JAY090029 - JAY090031</p>
 <p>JAY090032 - JAY090034</p>	 <p>JAY090035</p>	 <p>JAY090036</p>
 <p>JAY090037 - JAY090039</p>	 <p>JAY090040 - JAY090042</p>	 <p>JAY090043 - JAY090045</p>
 <p>JAY090046 - JAY090048</p>	 <p>JAY090049 - JAY090051</p>	 <p>JAY090052</p>
 <p>JAY090053</p>	 <p>JAY090054</p>	 <p>JAY090055</p>



JUP090047 & JUP090048



JCL020003




JCL020007




JCL030004



JCL030008



JCL030006



JCL030012



JCL030010



JUP030426 & JUP030417



JUP030427



JUP030428



JUP030429



JUP040201 & JUP040202



JUP040203 & JUP040204



JUP040205 & JUP040206



JUP040211 - JUP040214



JUP040207 - JUP040210



JUP110171



JUP110172



JUP090068



JUP110102



JUP110130



JUP110131



JUP110132



JUP110151



JUP110133



JUP110138



JUP110191 - JUP110198



JUP110181



V-Type




Knobby



JIV160101



JIV160103



JUP090067



JUP090061



JUP090505



JUP090506



JUP090507



JUP090508




JUP090509




JUP090510



JUP090511



JUP090040



JUP090069