

YOU-Q Wheelchairs

LUCA

QLASS

- Beautifully simple
- Exceptionally reliable
- Extremely manoeuvrable



John about Luca: "I like to get out and about to visit nice places. My Luca makes that possible and looks great as well: a chic one-piece chassis fitted with top-quality components. My Luca has rear-wheel drive and large wheels because I like riding around outdoors. But they can also supply the wheelchair with small wheels or mount the seating system the other way round on the chassis so it has front-wheel drive. That makes it even more manoeuvrable. The backrest and the seat are from the You-Q Qlass-collection. It's a fantastic seating system that can be perfectly adjusted to meet my requirements. In combination with the suspension, it provides me with optimal comfort. Most Lucas are grey; this one is as white as snow, but it has fun colour accents, which make it look really modern. It has given me back a large amount of the independence that I lost."

LUCA QCLASS



Chic design

The one-piece frame of the Luca is fitted with top-quality components. The Luca is very maintenance-friendly and has a chic, minimalist design. The colourful accents can be replaced in a trice.

Compact

The design of the Luca is well thought-out, which makes this wheelchair extremely compact and ideal for indoor use. Specially the front-wheel drive model is highly manoeuvrable.

Wide range

The Luca is available as an outdoor model with large wheels and an indoor model with smaller wheels. Both model types can be supplied in two speeds. You can also choose front-wheel or rear-wheel drive.

Perfect seating system

The Luca is fitted with the ultra-reliable Qclass seating system for optimal comfort.

DIMENSIONS		FWD	RWD
Maximum speed	km/h	6/10	6/10
Maximum user weight	kg	160	160
Batteries	Ah	50/60/78	50/60/78
Total length	mm	1100	1100
Total width			
12" drive wheels	mm	580	590
14" drive wheels	mm	610	625
Seat angle	°	0 - 7,5	0 - 7,5
Seat depth	mm	440 - 520	440 - 520
Seat width	mm	380 - 555	380 - 555
Seat height cushion excluded			
12" drive wheels	mm	425/450/475	425/450/475
14" drive wheels	mm	435/460/485	435/460/485
Backrest angle	°	89 - 118	89 - 118
Backrest height	mm	520 - 570	520 - 570
Diameter castor wheel			
8" castor wheels	mm	200	200
9" castor wheels	mm	225	225
Diameter drive wheel			
12" drive wheels	mm	320	320
14" drive wheels	mm	350	350
Powered backrest recline	°	88 - 120	88 - 120
Powered tilt in space Qclass	°	0 - 45	0 - 45
Powered seatlift ⁽¹⁾	mm	0 - 300	0 - 300
Powered legrests	°	6 - 70	6 - 70

1) adds 4.5 cm to the minimum seat height

ISO 7176-19
EN12182 (1999)
EN12184 (2009)

