

Mikra

Model No. **TE7A**

SHOPPRIDER



Special Features

- Easy to dismantle
- Twin detachable battery pack
- 12 amp power batteries
- Heaviest part 14.6 kg
- Tiller angle adjustment
- Swivel seat and flip up armrest
- Black non-marking puncture proof tyres
- Free wheel release handle
- Flip up adjustable armrests
- Also available in Blue



Specifications



| | | | |
|-----------------------------|---------------------------|--|------------------------------|
| Model No.* | TE7A | Battery Size 12v x 2 | 12 Ah |
| Class (Invalid Carriage) | 2 | Range** | 14-16 km / 8-10 miles |
| Overall Length | 1070 mm / 42" | Weight | 48.3kg / 107 lbs |
| Overall Width | 530 mm / 21" | Weight of Heaviest Part (Chassis Assembly) | 14.6kg / 32.2 lbs |
| Overall Height | 860 mm / 34" | Weight of Battery Pack | 6.5 kg (each) |
| Seat Width | 455 mm / 18" | Maximum Speed | 4 mph / 6.5 kmph |
| Seat Height (to ground) | 530-550mm / 21-22" | Wheel / Tyre Size | 200 mm / 8" |
| Armrest Height | 150 mm / 6" | Turning Radius | 1155 mm / 45.5" |
| Backrest Height (from Seat) | 350 mm / 14" | Battery Charger | 2 Amp Standard |
| Seat to Footrest / Floorpan | 368 mm / 14.5" | Ground Clearance*** | 89 mm / 3.5" |
| Safe Working Slope | 6° / 10% | | |
| Maximum User Weight | 113 kg / 18 st | | |

* Lower limb amputees should refer to user manual before use. ** Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel. *** Measurement to underside of chassis