

ROYALE 4 SPORT

Scooter

ROYALE

"The stylish canopy adds to this bold and dynamic design"



CE

The Royale range of scooters continues to evolve with the introduction of this state of the art Royale 4 Sport.

driveTM
www.drivemedical.co.uk

This stunning scooter features a stylish canopy, powerful 1500W motor and a maximum range up to an impressive 37 miles, providing the perfect balance of luxury comfort and performance.



Superb Comfort, Style & Performance...

- » Top speed of 8mph*.
 - » Maximum range of 28 miles (70Ah batts) or 37 miles (96Ah batts) on full battery charge*.
 - » Stylish modern design with canopy.
 - » Large 51cm (20") seat width accommodates larger users providing added comfort.
 - » Comfortable rotating seat with sliding and reclining facility.
 - » Height adjustable headrest to suit user requirements.
 - » Armrests are angle adjustable and flip up for easy side transfers.
 - » Automatic speed reduction when steering for added safety.
 - » All round suspension ensures a smooth ride and optimal handling.
 - » Tiller adjustment by lever to suit user requirements.
 - » Digital dash display with user friendly controls.
 - » High level indicators and lighting to improve outdoor visibility at dusk or when dark.
 - » Twin rear view mirrors for extra visibility.
 - » Large 36cm (14") front and 48cm (19") rear alloy wheels with black pneumatic tyres for a smoother ride and excellent ground clearance.
 - » Anti roll back system safety device prevents scooter from rolling backwards on hills.
 - » Freewheel facility allows movement of scooter without switching the motor on.
 - » Lap belt included for added safety.
 - » Highly powerful 1500W motor.
 - » Recommended MK batteries: 2 x 12V 79Ah or 2 x 12V 96Ah.
 - » Weight capacity 160kg (25st).
- * Subject to conditions of use.



Alloy wheels with black pneumatic tyres for a smoother ride

CODE: **HW006**

DESCRIPTION: **Royale 4 Sport**



Length:
169cm (67")



Height:
173cm (68")



Width:
75cm (30")



Max Capacity:
160kg (25st)



Ground Clearance:
12cm (5")



Max Gradient:
10°



High level indicators



Stylish design

Leaflet Code: **LL210**

If you would like to view this product please contact your local stockist, details below.